

ProCulture® Spinner Flasks

Designed to save time and increase cell production, ProCulture® Spinner Flasks have internal baffles for enhanced agitation and a magnetically driven vortex impeller that provides gentle top-to-bottom stirring.

- Impeller eliminates dead spots and extends above the media surface to increase aeration
- Sidearm openings large enough for 50 mL pipettes
- Flask base has central bump to prevent material from accumulating at the bottom center
- Entire assembly made of FDA approved materials
- Flask is Type 1, Class A borosilicate glass construction
- All components are autoclavable



ProCulture® Spinner Flasks

Complete Assembly

Catalog No.	Capacity	Center Neck Size	Side Arm Thread
BP-3400-100	125 mL	70mm	GL32
BP-3400-102	250 mL	70mm	GL32
BP-3400-104	500 mL	100mm	GL45
BP-3400-106	1 L	100mm	GL45
BP-3400-108	3 L	100mm	GL45
BP-3400-110	6 L	100mm	GL45
BP-3400-112	8 L	100mm	GL45
BP-3400-114	15 L	100mm	GL45
BP-3400-116	36 L	100mm	GL45



ProCulture® Wide-Mouth Spinner Flasks

Complete Assembly

- Provides easy component access assembly, cleaning, & sterilization

Catalog No.	Capacity	Center Neck Size	Side Arm Thread
BP-3500-108	3 L	120mm	GL45
BP-3500-110	6 L	120mm	GL45
BP-3500-112	8 L	120mm	GL45
BP-3500-114	15 L	120mm	GL45
BP-3500-116	36 L	120mm	GL45

ProCulture® Spinner Flasks with Shur-Loc® Caps

Complete Assembly

- Ideal for cultures that require a controlled environment
- Center neck cap has a unique compression O-ring that provides a visible seal

Catalog No.	Capacity	Center Neck Size	Side Arm Thread
BP-3600-108	3 L	140mm	GL45
BP-3600-110	6 L	140mm	GL45
BP-3600-112	8 L	140mm	GL45
BP-3600-114	15 L	140mm	GL45
BP-3600-116	36 L	140mm	GL45

ProCulture® is a registered trademark of



1172 NW Blvd., Vineland, NJ 08360 USA
1-800-220-5171 • Email: cs@wilmad-labglass.com
www.wilmad-labglass.com