

STOPPERS



Hollow Glass Stoppers w/ Std. Taper Joint

An Economy style standard taper stopper with an open bottom.
Made in the USA.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-10330-100	ST 14/20	1	\$22.81
LG-10330-102	ST 19/22	1	\$22.25
LG-10330-104	ST 24/40	1	\$34.73
LG-10330-105	ST 29/26	1	\$36.86
LG-10330-106	ST 29/42	1	\$37.28



Cylindrical Head Stoppers

Outer member stoppers.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-10370-100	ST 14/20	1	\$23.84
LG-10370-102	ST 19/22	1	\$25.83
LG-10370-104	ST 24/40	1	\$31.79
LG-10370-106	ST 29/42	1	\$43.31

STOPCOCKS



Balloon Adapters

These Air-Tite borosilicate glass adapters are designed with a T-bore stopcock for addition or venting inert atmospheres by attaching a balloon to one of two stems. Choose from different sized Standard Taper joints and plug type. IG-22 has a 4mm glass stopcock while IG-23T has a 4mm PTFE plug.

CATALOG NO.	DESCRIPTION	JOINT	QUANTITY / PACK	PRICE
IG-22-100	Glass Stopcock	ST 14/20	1	\$147.79
IG-22-104	Glass Stopcock	ST 24/40	1	\$159.45
IG-23T-100	PTFE Plug	ST 14/20	1	\$151.03
IG-23T-104	PTFE Plug	ST 24/40	1	\$157.89



Adapters with 90 Degree Bend to Hose Connection and Stopcock

Made of borosilicate glass, these adapters have a 90° bend to 5/16" O.D. hose connection on top and inner Standard Taper joint on bottom. Equipped with a stopcock, they are used to control gas/liquid addition in reaction systems. ML-210 and LG-1710 have glass plugs; ML-211T and LG-1710T have PTFE plugs.

CATALOG NO.	JOINT	PLUG TYPE	QUANTITY / PACK	PRICE
LG-1710-102	ST 24/40	Glass	1	\$73.04
LG-1710-104	ST 29/42	Glass	1	\$82.62
LG-1710T-202	ST 24/40	PTFE	1	\$99.46
LG-1710T-204	ST 29/42	PTFE	1	\$103.34
ML-210-700	ST 14/20	Glass	1	\$68.80
ML-210-702	ST 19/22	Glass	1	\$66.79
ML-211T-700	ST 14/20	PTFE	1	\$89.98
ML-211T-702	ST 19/22	PTFE	1	\$98.95

STOPCOCKS

Air-Tite Adapters with Vacuum Valve



These borosilicate glass adapters are equipped with vacuum valve for closing joints yet allowing for reconnection to other equipment with minimum exposure to air and vapors. The PTFE valve minimizes contamination and maintains a vacuum in the micron region. Standard Taper upper outer and bottom inner are same size.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
IG-202-100	ST 14/20	1	\$183.21
IG-202-102	ST 19/22	1	\$174.09
IG-202-104	ST 24/40	1	\$188.66

HOSE CONNECTORS

Air-Tite Hose Connector Adapter



This borosilicate glass adapter is supplied as a connection to LG-11071 and LG-11072 manifolds. Offered with size 15 O-ring joint for connection to heavy wall rubber tubing used with vacuum pumps and traps. Supplied with a LG-1022-530 Viton O-ring.

CATALOG NO.	APPROX OVERALL HEIGHT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-11073-100	50mm	3/8"	1	\$23.65



Straight Funnel Shaped Adapter

Made of borosilicate glass, this straight adapter has a funnel shaped interior that minimizes surface contact and two Standard Taper inner joints of the same size.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
ML-260-702	ST 19/22	1	\$84.99



Straight Adapters with Hose Connection and Inner Joint

With a Standard Taper inner joint at one end and a hose connection at the other, these straight adapters are made of borosilicate glass.

CATALOG NO.	HOSE CONNECTION O.D.	JOINT	QUANTITY / PACK	PRICE
LG-1461-102	3/8"	ST 24/40	1	\$28.67
LG-1461-104	3/8"	ST 29/42	1	\$32.02
ML-295-700	5/16"	ST 14/20	1	\$28.23
ML-295-702	5/16"	ST 19/22	1	\$29.53



Straight Adapters with Hose Connection and Ball Joint

With ball member of Socket Joint at one end and 5/16" hose connection at the other, these straight adapters are made of borosilicate glass.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1462-098	SJ 12/5	1	\$29.02
LG-1462-100	SJ 18/9	1	\$29.02
LG-1462-102	SJ 28/15	1	\$35.54
LG-1462-104	SJ 35/20	1	\$37.25
LG-1462-106	SJ 35/25	1	\$40.98



Straight Adapters with Hose Connection and Outer Joint

With a Standard Taper outer joint at one end and hose connection at other, these straight adapters are made of borosilicate glass.

CATALOG NO.	HOSE CONNECTION O.D.	JOINT	QUANTITY / PACK	PRICE
LG-1463-100	3/8"	ST 24/40	1	\$34.83
ML-296-700	5/16"	ST 14/20	1	\$34.80



Straight Adapters with Hose Connection and Socket Joint

With socket member of Ball and Socket Joint at one end and a 5/16" hose connection at the other, these straight adapters are made of borosilicate glass.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1464-098	SJ 12/5	1	\$31.62
LG-1464-100	SJ 18/9	1	\$37.02
LG-1464-102	SJ 28/15	1	\$37.69
LG-1464-104	SJ 35/25	1	\$45.55



90° Bend Adapters with Hose Connection and Ball Joint

These borosilicate glass adapters have a 90° bend to 5/16" hose connection on top and a ball member of Socket Joint at bottom.


CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1691-098	SJ 18/9	1	\$28.96
LG-1691-100	SJ 28/15	1	\$30.98
LG-1691-102	SJ 35/20	1	\$36.83
LG-1691-104	SJ 35/25	1	\$34.36



90° Bend Adapters with Hose Connection and Socket Joint

These borosilicate glass adapters have a 90° bend to 5/16" hose connection on top and a socket member of Ball and Socket Joint at other end.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1701-100	SJ 18/9	1	\$45.87
LG-1701-102	SJ 28/15	1	\$53.75
LG-1701-103	SJ 35/20	1	\$63.06



90° Bend Adapters with Hose Connection and Standard Taper Joint

These borosilicate glass adapters have a 90° bend to hose connector at top and a Standard Taper inner at bottom.

CATALOG NO.	HOSE CONNECTION O.D.	JOINT	QUANTITY / PACK	PRICE
LG-1702-104	3/8"	ST 29/42	1	\$29.56
ML-195-700	5/16"	ST 14/20	1	\$25.90
ML-195-702	5/16"	ST 19/22	1	\$27.90




90° Bend Full Length Adapters for Gas Inlet

These borosilicate glass gas inlet adapters have a 90° bend to 3/8" hose connector at top and a Standard Taper inner joint with a long stem sealed to its small end at bottom.

CATALOG NO.	JOINT	LENGTH BELOW JOINT	QUANTITY / PACK	PRICE
LG-1970-102	ST 24/40	200mm	1	\$91.14
LG-1970-104	ST 29/42	250mm	1	\$100.60

90° Bend Full Length Adapter with Drip Tip



These borosilicate glass adapters have a 90° bend to 3/8" hose connector at top and a Standard Taper inner with a short drip tip at bottom.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1980-100	ST 24/40	1	\$48.46

90° Air-Tite Adapters with Pressure/Vacuum Stopcock and Inner Joint



These Air-Tite borosilicate glass adapters have a 5/16" hose connector that is bent at a 90° angle on one end, and a Standard Taper inner at the other. A 2mm glass, pressure/vacuum stopcock lies between the two joints.

CATALOG NO.	JOINT	STOPCOCK BORE SIZE	QUANTITY / PACK	PRICE
IG-90-100	ST 14/20	2mm Glass	1	\$96.00
IG-90-104	ST 24/40	2mm Glass	1	\$105.52

105° Bend Adapter for Vacuum Distilling Systems



This borosilicate glass adapter is used in vacuum distillation systems and is bent at a 105° angle. Its top has a threaded Standard Taper 14/10 outer. The 8mm O.D. side hose connection is directly above the lower Standard Taper 14/10 - with drip tip extending below the joint. Includes 1 open top nylon screw cap and 1 Viton® O-ring.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
ML-390-600	ST 14/10	1	\$63.64

HOSE CONNECTORS



Serrated Hose Tube Connectors

Serrations are tooled in reverse on one end of tube of hose connectors so that the largest serration is at the end. The borosilicate glass tube is merely to act as a convenient means of holding the serrated portion while it is being sealed to an apparatus. When seal is completed, the tube is cut off and discarded. The diameter given is that of the smallest serration.

CATALOG NO.	DIAMETER	QUANTITY / PACK	PRICE
LG-10842-100	5/16"	1	\$4.80
LG-10842-102	3/8"	1	\$4.42
LG-10842-104	1/2"	1	\$4.76

HOSE CONNECTORS



Filtration Vacuum Adapters

Adapter for reduced pressure filtration using plain stem Buchner funnels. Top is tooled to accept special neoprene collar for funnels up to 4000mL capacity. Funnel sets in collar for vacuum seal and is easy to remove. See LG-1275C for collars.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1275-100	ST 24/40	3/8"	1	\$88.27
LG-1275-102	ST 29/42	3/8"	1	\$93.03
ML-265-700	ST 14/20	5/16"	1	\$82.62
ML-265-702	ST 19/22	5/16"	1	\$85.85

HOSE CONNECTORS



Adapter Collars for Vacuum Filtration Funnels

Neoprene adapter collars fit into reduced pressure filter adapters to accommodate plain stem Buchner funnels up to 3000mL in capacity. Funnels are set in collar for vacuum seal and are easy to remove. See LG-1275 and ML-265 for filter adapters.

CATALOG NO.	FITS FUNNEL SIZE	QUANTITY / PACK	PRICE
LG-1275C-100	15 and 60mL	1	\$8.61
LG-1275C-102	600mL	1	\$15.66
LG-1275C-104	3000mL	1	\$16.45



Fritted Adapters with Hose Connection

Made of borosilicate glass, these adapters have an extra coarse disc between top hose connection and bottom inner joint to prevent introduction of solids into a reaction mixture.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-7152-100	ST 14/20	1	\$49.46
LG-7152-102	ST 19/22	1	\$46.37
LG-7152-104	ST 24/40	1	\$52.39
LG-7152-108	ST 29/42	1	\$53.39



Inert Gas Adapter

This borosilicate glass adapter has an inner Standard Taper 24/40 at bottom and two 90 degree hose connections at top that allow inert gas flow over apparatus to maintain an oxygen free system.

CATALOG NO.	HOSE CONNECTION O.D.	JOINT	QUANTITY / PACK	PRICE
LG-1465-100	3/8"	ST 24/40	1	\$39.12



Hollow Glass Stoppers

An Economy style interchangeable taper-ground stopper. Made in the USA.

CATALOG NO.	STOPPER NUMBER	QUANTITY / PACK	PRICE
LG-10326-108	22	1	\$22.88
LG-10326-110	27	1	\$27.45



Glass and PTFE Straight Adapters with Hose Connection

These borosilicate glass adapters feature a 5/16" O.D. hose connection at top and Standard Taper 24/40 on bottom. With a stopcock between joints, model LG-1711 has a 2mm glass plug and LG-1711T has a 2mm PTFE plug.

CATALOG NO.	JOINT	STOPCOCK BORE SIZE	QUANTITY / PACK	PRICE
LG-1711-102	ST 24/40	2mm Glass	1	\$81.32
LG-1711T-202	ST 24/40	2mm PTFE	1	\$100.02



Fritted Adapter with PTFE Stopcock

Supplied with a coarse porosity disc, this borosilicate glass adapter features a 5/16" O.D. hose connection at top and Standard Taper 24/40 on bottom. The stopcock between the joints has a PTFE plug with a 2mm bore.

CATALOG NO.	JOINT	STOPCOCK BORE SIZE	QUANTITY / PACK	PRICE
LG-7151T-104	ST 24/40	2mm PTFE	1	\$126.45

HOSE CONNECTORS



Straight Vacuum Adapters

For vacuum filtration into vessels. 3/8" hose connector is provided for attachment to vacuum source. Top opening accommodates a #6 stopper.

CATALOG NO.	JOINT	STOPPER NUMBER	QUANTITY / PACK	PRICE
LG-1272-100	ST 24/40	6	1	\$70.96
LG-1272-102	ST 29/42	6	1	\$74.81



Straight Vacuum Adapters with 3/8 Inch Hose Connection

With a standard taper outer joint at one end and standard taper inner joint at the other; hose connection for vacuum source.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1330-100	ST 24/40 Inner, ST 24/40 Outer	3/8"	1	\$60.87
LG-1330-102	ST 29/42 Inner, ST 29/42 Outer	3/8"	1	\$67.09

HOSE CONNECTORS

Bent and Straight Hose Connections

Borosilicate glass hose connecting units are used to connect to apparatus with GL14 or GL18 thread sizes. Complete connecting unit consists of connector - either straight or bent - and cap with hole in center.



CATALOG NO.	DESCRIPTION	GL THREAD SIZE	QUANTITY / PACK	PRICE
LG-10818-100	Complete Unit - Bent	GL14	1	\$4.69
LG-10818-102	Complete Unit - Straight	GL14	1	\$4.94
LG-10818-104	Complete Unit - Bent	GL18	1	\$12.13
LG-10818-106	Complete Unit - Straight	GL18	1	\$12.73
LG-10818-120	Bent Connector	GL14	1	\$2.54
LG-10818-122	Straight Connector	GL14	1	\$2.38
LG-10818-124	Bent Connector	GL18	1	\$10.42
LG-10818-126	Straight Connector	GL18	1	\$10.13

T-Shaped Tubing Connectors

Borosilicate glass tube connectors are T-shaped with serrated ends for connection to tubing.



CATALOG NO.	O.D.	QUANTITY / PACK	PRICE
LG-10825-100	8mm	1	\$20.90
LG-10825-102	10mm	1	\$23.13


HOSE CONNECTORS

105 Degree Bent Vacuum Adapter

Adapter with a drip tip and a 3/8" hose connection.



CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1371-100	ST 24/40	1	\$103.73



Straight Vacuum Adapters with Hose Connection and Sealed Delivery Tube

With a Standard Taper outer joint at top and a Standard Taper inner joint of the same size at bottom, 3/8" side hose connection and sealed delivery tube.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1840-100	ST 24/40 Outer, ST 24/40 Inner	3/8"	1	\$74.35
LG-1840-102	ST 29/42 Outer, ST 29/42 Inner	3/8"	1	\$80.43
ML-375-700	ST 14/20 Outer, ST 14/20 Inner	3/8"	1	\$58.28
ML-375-702	ST 19/22 Outer, ST 19/22 Inner	3/8"	1	\$54.10

HOSE CONNECTORS



Y-Shaped Tubing Connectors

Borosilicate glass tube connectors are Y-shaped with serrated ends for connection to tubing.

CATALOG NO.	O.D.	QUANTITY / PACK	PRICE
LG-10835-100	8mm	1	\$25.17
LG-10835-102	10mm	1	\$25.65


HOSE CONNECTORS



Straight Vacuum Adapters with 5/16 Inch Hose Connection

Sealed in delivery tube extends 90mm below lower inner joint. Upper and lower joints are the same. Supplied with 5/16" hose connection.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1842-100	ST 24/40	5/16"	1	\$82.04
LG-1842-102	ST 29/42	5/16"	1	\$78.83
ML-380-700	ST 14/20	5/16"	1	\$58.04
ML-380-702	ST 19/22	5/16"	1	\$63.46



90 Degree Bent Gas Inlet Vacuum Adapters

Provided with a standard taper inner joint and delivery tube extending 250mm below lower joint. 3/8" hose connectors allow attachment of flexible tubing.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-2050-102	ST 24/40	3/8"	1	\$108.34
LG-2050-104	ST 29/42	3/8"	1	\$116.07

HOSE CONNECTORS



Air-Tite Adapters with Pressure/Vacuum Stopcock and Outer Joint

Made of borosilicate glass, these Air-Tite adapters are used for evacuation and addition of inert gases. They consist of a Standard Taper outer joint on bottom, 2mm glass stopcock, and 5/16" hose connector on top.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
IG-10-100	ST 14/20	1	\$128.66
IG-10-104	ST 24/40	1	\$113.58

STRAIGHT ADAPTERS



Anti-Splash Adapters

Made of borosilicate glass, these anti-splash adapters are equipped with Standard Taper joints at top (outer) and bottom (inner). When adapter is inserted into the neck of a flask the guard prevents splash up. Overall height is minimized to provide best operation, especially during distillation.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1242-098	ST 14/20	1	\$65.98
LG-1242-102	ST 24/40	1	\$78.95



Anti-Splash Adapters with Coarse Frit

Made of borosilicate glass, these anti-splash adapters with a coarse fritted glass are equipped with Standard Taper joints at top (outer) and bottom (inner). When adapter is inserted into the neck of a flask the guard prevents splash up. Overall height is minimized to provide best operation, especially during distillation.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1243-098	ST 14/20	1	\$82.83
LG-1243-100	ST 24/25	1	\$89.19
LG-1243-102	ST 24/40	1	\$94.34

Bushing Style Reducing Adapters

This reducing adapter provides compact Standard Taper assemblies. The larger Standard Taper size is on the outside bottom of bushing while the smaller Standard Taper dimension is on top within the bore of bushing.



CATALOG NO.	TOP JOINT	BOTTOM JOINT	QUANTITY / PACK	PRICE
LG-1240-106	ST 10/30	ST 24/40	1	\$50.10
LG-1240-108	ST 14/35	ST 24/40	1	\$46.11
LG-1240-110	ST 19/38	ST 24/40	1	\$43.92
LG-1240-112	ST 10/30	ST 29/42	1	\$52.47
LG-1240-114	ST 14/35	ST 29/42	1	\$54.04
LG-1240-116	ST 19/38	ST 29/42	1	\$48.64
LG-1240-118	ST 24/40	ST 29/42	1	\$49.82
LG-1240-120	ST 24/40	ST 34/45	1	\$65.89
LG-1240-122	ST 29/42	ST 34/45	1	\$82.53
LG-1240-126	ST 24/40	ST 45/50	1	\$110.94
LG-1240-128	ST 29/42	ST 45/50	1	\$104.81
LG-1240-130	ST 34/45	ST 45/50	1	\$108.63
LG-1240-133	ST 29/42	ST 55/50	1	\$138.06
LG-1240-136	ST 45/50	ST 55/50	1	\$136.54
ML-60-700	ST 7/12	ST 14/20	1	\$35.77
ML-60-702	ST 10/18	ST 14/20	1	\$33.59
ML-60-712	ST 10/30	ST 19/22	1	\$43.17
ML-60-714	ST 14/20	ST 19/22	1	\$38.20

Reducing and Enlarging Adapters



These borosilicate glass reducing and enlarging adapters allow connection of dissimilar Standard Taper joints. Outer joint is at the top and inner joint is at the bottom.

CATALOG NO.	TOP JOINT	BOTTOM JOINT	QUANTITY / PACK	PRICE
LG-1190-112	ST 10/30	ST 19/38	1	\$37.40
LG-1190-114	ST 10/30	ST 24/40	1	\$42.91
LG-1190-162	ST 19/38	ST 24/40	1	\$45.37
LG-1190-178	ST 24/40	ST 24/40	1	\$52.94
LG-1190-180	ST 24/40	ST 29/42	1	\$47.30
LG-1190-182	ST 24/40	ST 34/45	1	\$54.38
LG-1190-184	ST 24/40	ST 40/50	1	\$76.58
LG-1190-186	ST 24/40	ST 45/50	1	\$71.96
LG-1190-190	ST 24/40	ST 55/50	1	\$95.93
LG-1190-198	ST 29/42	ST 34/45	1	\$68.56
LG-1190-206	ST 29/42	ST 55/50	1	\$103.12
LG-1200-100	ST 19/38	ST 10/30	1	\$40.59
LG-1200-106	ST 24/40	ST 10/30	1	\$42.18
LG-1200-110	ST 24/40	ST 14/35	1	\$43.75
LG-1200-112	ST 24/40	ST 19/38	1	\$45.12
LG-1200-118	ST 29/42	ST 19/38	1	\$50.32
LG-1200-120	ST 29/42	ST 24/40	1	\$50.31
LG-1200-124	ST 34/45	ST 24/40	1	\$57.88
LG-1200-126	ST 34/45	ST 29/42	1	\$61.45
LG-1200-136	ST 45/50	ST 24/40	1	\$102.60
LG-1200-140	ST 55/50	ST 24/40	1	\$108.14
ML-235-750	ST 29/26	ST 19/22	1	\$48.11
ML-235-706	ST 14/20	ST 19/22	1	\$35.51
ML-235-712	ST 14/20	ST 24/40	1	\$40.23
ML-235-720	ST 19/22	ST 14/20	1	\$38.69
ML-235-726	ST 19/22	ST 24/40	1	\$42.69
ML-235-730	ST 19/22	ST 29/42	1	\$44.78
ML-235-744	ST 24/40	ST 14/20	1	\$41.09
ML-235-746	ST 24/40	ST 19/22	1	\$41.92



Straight Adapters with Standard Taper Joints

These borosilicate glass straight adapters have same size inner Standard Taper joint at each end.

CATALOG NO.	LENGTH BETWEEN JOINTS	JOINT	QUANTITY / PACK	PRICE
LG-1651-100	30mm	ST 24/40	1	\$58.77
LG-1651-102	75mm	ST 24/40	1	\$53.43
LG-1651-104	120mm	ST 24/40	1	\$56.35
LG-1651-106	170mm	ST 24/40	1	\$59.86
LG-1651-202	75mm	ST 29/42	1	\$56.08
LG-1651-204	120mm	ST 29/42	1	\$56.08
LG-1651-206	170mm	ST 29/42	1	\$57.39



Straight Connecting Adapters

These borosilicate glass straight adapters feature a Standard Taper outer joint at one end and the same size Standard Taper inner joint at the other.

CATALOG NO.	JOINT	LENGTH BETWEEN JOINTS	QUANTITY / PACK	PRICE
LG-1290-100	ST 24/40	30mm	1	\$48.75
LG-1290-102	ST 24/40	75mm	1	\$48.90
LG-1290-104	ST 24/40	120mm	1	\$55.94
LG-1290-106	ST 24/40	170mm	1	\$57.03
LG-1290-110	ST 29/42	30mm	1	\$48.90
LG-1290-112	ST 29/42	75mm	1	\$52.47
LG-1290-114	ST 29/42	120mm	1	\$59.14
ML-245-700	ST 14/20	50mm	1	\$36.51

STRAIGHT ADAPTERS



Open Top Bushing/Threaded Glass Connectors

These borosilicate glass connectors are open-ended, have external threads, and are available in Nylon and PTFE. They are for use with threaded glass connectors (LG-2075) for connection to reducing end tubes. One Viton O-ring included.

CATALOG NO.	ACCEPTS TUBING O.D.	THREADED CONNECTOR SIZE	MATERIAL TYPE	QUANTITY / PACK	PRICE
LG-2079-100	6-7mm	GL 7	Nylon	1	\$3.92
LG-2079-102	9-10.5mm	GL 11	Nylon	1	\$30.52
LG-2079-104	12.5-14mm	GL 15	Nylon	1	\$23.01
LG-2079-200	6-7mm	GL 7	PTFE	1	\$17.91
LG-2079-202	9-10.5mm	GL 11	PTFE	1	\$45.42
LG-2079-204	12.5-14mm	GL 15	PTFE	1	\$50.88

STRAIGHT ADAPTERS



Flexible PTFE Adapters with Standard Taper Joint

These flexible PTFE connectors relieve stress and allow for movement of glass joint connection and mis-alignment. Top and bottom Standard Taper joint connections are same size.

CATALOG NO.	HEIGHT	JOINT	QUANTITY / PACK	PRICE
LG-1014T-200	115mm	ST 14/20	1	\$76.34
LG-1014T-202	90mm	ST 19/22	1	\$79.98
LG-1014T-204	120mm	ST 19/38	1	\$82.89
LG-1014T-206	130mm	ST 24/40	1	\$103.53
LG-1014T-208	140mm	ST 29/42	1	\$108.77
LG-1014T-210	145mm	ST 34/45	1	\$139.40
LG-1014T-212	155mm	ST 45/50	1	\$192.58

BENT ADAPTERS



75° Bend Adapters

These borosilicate glass adapters are bent at a 75° angle and feature a Standard Taper inner joint of the same size at each end.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1390-100	ST 24/40	1	\$47.93
LG-1390-102	ST 29/42	1	\$49.26
ML-20-700	ST 14/20	1	\$33.31
ML-20-702	ST 19/22	1	\$33.05



Impinger 90° Adapters

Bent at a 90° angle, these borosilicate glass adapters feature a Socket joint at one end and a Ball joint of the same size at the other.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-6889A-100	SJ 12/5	1	\$42.04
LG-6889A-102	SJ 18/9	1	\$44.60
LG-6889A-104	SJ 28/15	1	\$49.27

105° Bend Adapters



These borosilicate glass adapters are bent at a 105° angle and feature a Standard Taper inner joint on one end and a Standard Taper outer joint of the same size at other end.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1370-100	ST 24/40	1	\$49.27
LG-1370-102	ST 29/42	1	\$54.04
ML-35-700	ST 14/20	1	\$49.84
ML-35-702	ST 19/22	1	\$54.31

105° Bend Adapters with Drip Tip



These borosilicate glass adapters have an outer Standard Taper joint at one end and a small tube at the other, angled to form a drip tip. Bend gives an included angle of 105°.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1350-100	ST 24/40	1	\$38.83
LG-1350-102	ST 29/42	1	\$38.66
ML-10-700	ST 14/20	1	\$34.42
ML-10-702	ST 19/22	1	\$32.25

105° Bend Adapters for Distilling Heads



These borosilicate glass bent adapters are used to connect distilling head, LG-6171 to distilling receiver, LG-6381. Bent at 105° angle, adapters have a Socket Joint to connect takeoff tube on the distilling head and the lower portion has a Ball Joint for attachment to receiver.


CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-6190-100	SJ 18/9	1	\$52.59
LG-6190-102	SJ 28/15	1	\$55.91



Vacuum Adapters with 105 Degree Bend, Style B, Standard Taper Joints

These borosilicate glass vacuum-type adapters feature Standard Taper joints (top outer and bottom inner) and a sealed-in inner delivery tube on bottom joint which ends in a slanted drip tip. Designed with a side hose connector (Style A is pointed away from adapter) for connection to a vacuum source.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
ML-390-700	ST 14/20 top outer, bottom inner	1	\$62.20
ML-390-702	ST 19/22 top outer, bottom inner	1	\$61.76



105° Bend Adapters with Drip Tip and Inner/Outer Joint

Made of borosilicate glass, these vacuum adapters feature a delivery tube with a drip tip that extends below the lower inner Standard Taper joint. A Standard Taper outer joint of the same size is at other end of the 105° bent adapter. Connection to vacuum source is provided by a hose connector.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1860-102	ST 24/40	1	\$72.22
LG-1860-104	ST 29/42	1	\$75.52


Connecting Adapters

These borosilicate glass connecting tube adapters have two parallel Standard Taper inner joints that connect at a 75° angle to an upper outer joint.



CATALOG NO.	OFFSET CENTER TO CENTER	JOINT	QUANTITY / PACK	PRICE
LG-1760-100	170mm	ST 24/40	1	\$62.61
LG-1760-102	200mm	ST 29/42	1	\$73.83

DISTILLATION ADAPTERS



Distilling Adapters with 75 Degree Angle Side Arm

With sidearm at angle of 75° from lower part of main tube, these borosilicate glass distilling adapters have three Standard Taper joints of the same size. Sidearm and lower tube joint are inner and upper joint is outer.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1470-100	ST 24/40	1	\$71.21
LG-1470-102	ST 29/42	1	\$142.19

DISTILLATION ADAPTERS



Bump Adapter Evaporator Trap

With an extra coarse frit sealed between the outer upper joint and the lower inner joint this adapter prevents the contents in the flask from entering into the condenser when "bumping" occurs. All traps have a 24/40 top outer joint.


CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-6877-100	ST 14/20 Lower	1	\$98.82
LG-6877-102	ST 24/40 Lower	1	\$108.67
LG-6877-106	ST 29/42 Lower	1	\$114.76



Bump Pear Shaped RotoVap Trap

Improved design returns condensate to evaporator flask. Pear shaped flask slope is similar to rotary evaporator angle.

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-6883-100	250mL	ST 24/40 Upper, ST 14/20 Lower	1	\$93.10
LG-6883-102	250mL	ST 24/40 Upper, ST 19/22 Lower	1	\$95.06
LG-6883-104	250mL	ST 24/40 Upper, ST 24/40 Lower	1	\$94.04



Self Washing Evaporator Trap

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-6878-100	100mL	ST 24/40 to 14/20	1	\$133.83



Rotary Evaporator Traps

The upper part of trap has a standard taper outer joint.

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-6880-100	100mL	ST 24/40 Upper, ST 14/20 Lower	1	\$135.89
LG-6880-102	100mL	ST 24/40 Upper, ST 19/22 Lower	1	\$134.44
LG-6880-104	100mL	ST 24/40 Upper, ST 24/40 Lower	1	\$134.95
LG-6880-105	100mL	ST 29/26 Upper, ST 29/26 Lower	1	\$148.53
LG-6880-106	250mL	ST 24/40 Upper, ST 14/20 Lower	1	\$139.47
LG-6880-108	250mL	ST 24/40 Upper, ST 19/22 Lower	1	\$142.48
LG-6880-110	250mL	ST 24/40 Upper, ST 24/40 Lower	1	\$138.45
LG-6880-114	250mL	ST 29/26 Upper, ST 29/26 Lower	1	\$142.70
LG-6880-115	250mL	ST 24/40 Upper, ST 29/26 Lower	1	\$143.10
LG-6880-116	500mL	ST 24/40 Upper, ST 24/40 Lower	1	\$154.54
LG-6880-117	500mL	ST 24/40 Upper, ST 29/26 Lower	1	\$145.47
LG-6880-118	500mL	ST 29/26 Upper, ST 29/26 Lower	1	\$147.70
LG-6880-120	1000mL	ST 24/40 Upper, ST 29/42 Lower	1	\$157.71



Modified Rotary Evaporator Traps

Similar to LG-6880 except inner vapor tube is sealed to top of trap to prevent entry into the condenser. All top joints are ST 24/40.

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-6880A-100	100mL	ST 14/20 Lower	1	\$141.66
LG-6880A-102	100mL	ST 24/40 Lower	1	\$132.90
LG-6880A-104	250mL	ST 14/20 Lower	1	\$133.52
LG-6880A-106	250mL	ST 24/40 Lower	1	\$137.94
LG-6880A-108	250mL	ST 29/42 Lower	1	\$138.74
LG-6880A-110	500mL	ST 14/20 Lower	1	\$134.15
LG-6880A-112	500mL	ST 24/40 Lower	1	\$137.94
LG-6880A-114	500mL	ST 29/42 Lower	1	\$138.74



Rotary Evaporator Traps with Two 3mm Holes

Base of inner tube has two 3mm openings to allow drain-back of foam condensate to the evaporator Flask. Top joint is a Standard Taper outer.

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-6881-100	100mL	ST 24/40 Upper, ST 14/20 Lower	1	\$142.48
LG-6881-102	100mL	ST 24/40 Upper, ST 19/22 Lower	1	\$152.13
LG-6881-104	100mL	ST 24/40 Upper, ST 24/40 Lower	1	\$143.48
LG-6881-105	100mL	ST 29/26 Upper, ST 29/26 Lower	1	\$141.51
LG-6881-106	250mL	ST 24/40 Upper, ST 14/20 Lower	1	\$135.78
LG-6881-108	250mL	ST 24/40 Upper, ST 19/22 Lower	1	\$154.47
LG-6881-110	250mL	ST 24/40 Upper, ST 24/40 Lower	1	\$147.12
LG-6881-114	500mL	ST 24/40 Upper, ST 19/22 Lower	1	\$157.08
LG-6881-115	250mL	ST 24/40 Upper, ST 29/26 Lower	1	\$157.06
LG-6881-116	500mL	ST 24/40 Upper, ST 24/40 Lower	1	\$157.06
LG-6881-117	500mL	ST 24/40 Upper, ST 29/26 Lower	1	\$149.49
LG-6881-118	1000mL	ST 24/40 Upper, ST 24/40 Lower	1	\$151.18
LG-6881-120	1000mL	ST 24/40 Upper, ST 29/42 Lower	1	\$151.68



Rotary Evaporator Connectors and Adapters

Adapter has a top outer Standard Taper joint of the listed size and a 20-400 G.P.I. thread on the bottom. A PTFE lined, polypropylene double threaded connector attaches a standard laboratory vial.

CATALOG NO.	DESCRIPTION	JOINT	THREADED CONNECTOR SIZE	QUANTITY / PACK	PRICE
LG-6879-110	Connector		13-425	1	\$19.73
LG-6879-112	Connector		20-400	1	\$24.86
LG-6879-116	Connector		24-400	1	\$35.34
LG-6879-120	Glass Adapter Only	ST 24/40		1	\$55.52
LG-6879-122	Glass Adapter Only	ST 29/26		1	\$78.24

INLET ADAPTERS



Universal Inlet Adapters with Hose Connector, Open Screw Cap

These borosilicate glass adapters consist of a hose connector for attachment to vacuum/gas source (diameter accommodation range of 2 to 17mm), top open compression cap (fits thermometers, gas dispersion tubes, pH probes) to seal inlet using Viton® O-ring (included), and bottom inner Standard Taper 24/40. NOTE: Top opening can also be sealed using PTFE faced silicone septum (sold separately).

CATALOG NO.	ACCOMODATION RANGE	O-RING CODE	QUANTITY / PACK	PRICE
LG-1452-300	2-5mm	LG-1022-527	1	\$53.21
LG-1452-304	5-8mm	LG-1022-519	1	\$53.21
LG-1452-306	8-11mm	LG-1022-522	1	\$79.47



Straight Thermometer Adapters

With Standard Taper inner member at one end and plain tubulature at opposite, these borosilicate glass straight adapters are used for attachment of a thermometer, gas inlet, bleed tube or stirrer. A rubber fitment (LG-1988) is supplied to connect these pieces to the adapter.

CATALOG NO.	DESCRIPTION	BOTTOM JOINT	QUANTITY / PACK	PRICE
ML-280-700	Adapter w/ Rubber Fitment	ST 14/20	1	\$23.79
ML-280-702	Adapter w/ Rubber Fitment	ST 19/22	1	\$22.77
LG-1985-102	Adapter w/ Rubber Fitment	ST 24/40	1	\$36.10
LG-1985-103	Adapter w/ Rubber Fitment	ST 29/26	1	\$38.38
LG-1985-104	Adapter w/ Rubber Fitment	ST 29/42	1	\$42.18
LG-1988-100	Replacement Rubber Fitment		1	\$3.61



Threaded Seal Thermometer Adapters

Made of borosilicate glass, these straight adapters feature ground joint at bottom and 7mm GL14 threaded compression fitting at top. They are for use with thermometers, bleed tubes, inlet tubes, etc. when used with threaded cap (LG-10818) and PTFE faced silicone seal ring (LG-14038) – each sold separately.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1456-100	ST 24/40	1	\$36.56
LG-1456-104	ST 29/42	1	\$36.68



Off-Centered Thermometer Adapters

These borosilicate glass off-centered adapters feature an inner Standard Taper joint at one end and an offset Standard Taper 10/30 outer joint at the other. Offset joint is set on approximate 10° angle to inner. Thermometer bulb position within apparatus may be changed by rotating adapter.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1983-102	ST 24/40	1	\$57.02
LG-1983-104	ST 29/42	1	\$61.63



Three-Way Thermometer Adapters, 75 Degree Angle Side Arm

These borosilicate glass three-way thermometer adapters have a 7mm threaded fitting at top, a bottom inner Standard Taper joint, and 75° angle inner side joint of same size. Supplied with cap and O-ring for inert and greaseless thermometer connection and easy adjustment of thermometer position.

CATALOG NO.	ACCOMODATION RANGE	APPROX OVERALL HEIGHT	APPROX OVERALL WIDTH	JOINT	QUANTITY / PACK	PRICE
LG-1521-098	5-8mm	110mm	55mm	ST 19/22	1	\$80.26
LG-1521-102	5-8mm	165mm	110mm	ST 29/42	1	\$163.34



Universal Inlet Adapter with Hose Connector, Open Top Nylon Screw Cap

Made of borosilicate glass, this Universal adapter has Standard Taper 14/10 inner joint on bottom and 13-425 screw thread on top that forms an O-ring compression seal with thermometers, bleed tubes, inlet tubes and accommodates 5.0 - 7.5mm O.D. tubes or thermometers. Side hose connection is 8mm O.D. Includes one 20-400 open top nylon screw cap and one size 014 Viton® O-ring.

CATALOG NO.	QUANTITY / PACK	PRICE
LG-1452-090	1	\$48.65



Distilling Adapter with 75 Degree Angle Side Arm

This borosilicate glass three-way adapter has a Standard Taper 10/30 outer joint on top to accommodate a thermometer, bottom inner Standard Taper 24/40, and 75° angle inner side joint of same size.

CATALOG NO.	JOINT	IMMERSION	THERMOMETER JOINT	QUANTITY / PACK	PRICE
LG-1530-100	ST 24/40	75mm	ST 10/30	1	\$117.18



Distilling Adapters with 75 Degree Angle Extended Side Arm

These borosilicate glass three-way adapters have a Standard Taper outer joint on top to accommodate a thermometer, bottom Standard Taper inner, and 75° angled sidearm with Standard Taper inner. Angled sidearm is extended below main tube and bent straight down so both bottom and side joints are parallel and of same size.

CATALOG NO.	JOINT	IMMERSION	THERMOMETER JOINT	OFFSET CENTER TO CENTER	QUANTITY / PACK	PRICE
LG-1740-100	ST 24/40	75mm	ST 10/30	160mm	1	\$101.57
LG-1740-102	ST 29/42	75mm	ST 10/30	175mm	1	\$106.58

Universal Inlet Adapters with Screw Cap



Borosilicate glass universal inlet adapters are capable of sealing thermometers, syringes, pipets, etc. Insert probes can be positioned within the apparatus by loosening the screw cap on top of adapter and tightening to create a seal. Bottom of adapter features a Standard Taper joint. Supplied as a complete assembly with a screw cap and Viton® O-ring.

CATALOG NO.	JOINT	ACCOMODATION RANGE	QUANTITY / PACK	PRICE
LG-1451-100	ST 14/10	2-5mm	1	\$39.52
LG-1451-090	ST 14/10	5-5.7mm	1	\$33.87
LG-1451-102	ST 14/10	5.5-6.5mm	1	\$32.71
LG-1451-104	ST 14/20	2-5mm	1	\$57.81
LG-1451-108	ST 14/20	5-8mm	1	\$41.82
LG-1451-200	ST 19/22	2-5mm	1	\$36.60
LG-1451-204	ST 19/22	5-8mm	1	\$42.73
LG-1451-300	ST 24/40	2-5mm	1	\$40.15
LG-1451-304	ST 24/40	5-8mm	1	\$42.49
LG-1451-306	ST 24/40	8-11mm	1	\$50.55
LG-1451-308	ST 24/40	11-14mm	1	\$54.27
LG-1451-400	ST 29/42	2-5mm	1	\$41.82
LG-1451-404	ST 29/42	5-8mm	1	\$44.24
LG-1451-406	ST 29/42	8-11mm	1	\$58.47

Joint Converting Thermometer Adapters



With a standard taper joint at the bottom and top this adapter fits thermometers calibrated for 75mm immersion. The bulb will extend 40mm below the bottom of the adapter.

CATALOG NO.	TOP JOINT	BOTTOM JOINT	QUANTITY / PACK	PRICE
ML-290-702	ST 10/18	ST 14/20	1	\$48.27
ML-290-704	ST 7/12	ST 19/22	1	\$51.05
ML-290-706	ST 10/18	ST 19/22	1	\$49.01

KJELDAHL ADAPTERS



Kjeldahl Adapter with Delivery Tube

This borosilicate glass Kjeldahl adapter has a Standard Taper inner joint with drip tip at one end just below the round trap. A 150mm long thick-walled delivery tube angled at 75° is at other end.

CATALOG NO.	JOINT	TUBING O.D.	QUANTITY / PACK	PRICE
LG-1660-100	ST 24/40	8mm	1	\$81.51



Kjeldahl Adapters with Trap

This Kjeldahl trap borosilicate glass adapter has two parallel Standard Taper inner joints of the same size approximately 200mm apart, center-to-center. Connecting tube is at a 75° angle from joint with drip tip below round trap.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1662-100	ST 24/40	1	\$183.63
LG-1662-102	ST 29/26	1	\$183.63
LG-1662-104	ST 29/42	1	\$183.63



Kjeldahl Distilling Trap Straight Adapters

With upper Standard Taper outer and same size lower Standard Taper inner, these borosilicate glass Kjeldahl adapters are placed above flask to prevent foaming-type reactions from entering either distillation column or condenser. Approximate bulb capacity is 125mL.

CATALOG NO.	TOP JOINT	BOTTOM JOINT	QUANTITY / PACK	PRICE
LG-1910-100	ST 24/40	ST 24/40	1	\$64.96
LG-1910-102	ST 29/42	ST 29/42	1	\$81.06
LG-1910-104	ST 24/40	ST 14/20	1	\$75.07
LG-1910-106	ST 24/40	ST 19/22	1	\$74.27
LG-1910-108	ST 29/42	ST 14/20	1	\$72.99
LG-1910-110	ST 29/42	ST 19/22	1	\$75.14

OFFSET ADAPTERS



Off-Set Adapters

These borosilicate glass off-set adapters have a Standard Taper outer joint on one end and same size Standard Taper inner joint on the other.

CATALOG NO.	JOINT	OFFSET CENTER TO CENTER	QUANTITY / PACK	PRICE
LG-1580-100	ST 24/40	22mm	1	\$51.23
LG-1580-102	ST 29/42	28mm	1	\$77.39

MULTI-WAY ADAPTERS



Full Length Three-Way Adapters

With sidearm at angle of 105° from lower tube, these borosilicate glass three-way adapters have same size Standard Taper joints - outer on sidearm and upper tube, inner on lower tube.


CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1440-096	ST 14/20	1	\$75.02
LG-1440-100	ST 24/40	1	\$79.42
LG-1440-102	ST 29/42	1	\$121.81



Three-Way Adapter with Ball and Socket Joint

This borosilicate glass three-way adapter has Socket Joints of the same size at ends.


CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1450-100	SJ 28/15	1	\$43.46



Claisen Adapters

With one vertical sidearm and one at 45° angle to the main stem, these borosilicate glass four-way adapters have Standard Taper joints of the same size. Claisen adapters allow an additional entry into single or multiple neck flasks and are useful for reflux with addition reactions.


CATALOG NO.	JOINT	OFFSET CENTER TO CENTER	QUANTITY / PACK	PRICE
LG-1812-100	ST 24/40	100mm	1	\$142.79
LG-1812-102	ST 29/42	100mm	1	\$185.41



PTFE Thermometer Adapters with Compression Cap

These PTFE inner joint adapters have a threaded compression cap, four saddle O-rings and one Viton® O-ring for use with plain stem thermometers and inlet tubes. Accepts thermometers with O.D.'s from 3 to 7mm.

CATALOG NO.	JOINT	FITS OUTER JOINTS	QUANTITY / PACK	PRICE
LG-1453-100	ST 10/18	ST 10/18, ST 10/30	1	\$25.23
LG-1453-102	ST 14/20	ST 14/20, ST 14/35	1	\$30.29
LG-1453-104	ST 19/22	ST 19/22, ST 19/38	1	\$31.78
LG-1453-106	ST 24/25	ST 24/25, ST 24/40	1	\$37.30
LG-1453-108	ST 29/26	ST 29/26, ST 29/42	1	\$42.73
LG-1453-110	ST 34/28	ST 34/28, ST 34/45	1	\$61.77



Claisen Adapter with Caps

Made of borosilicate glass, this three-way adapter consists of upper outer joints with screw caps and a lower inner joint. All joints are Standard Taper 14/10. Includes two open top nylon screw caps, two Viton® O-rings, and two PTFE-faced silicone septa. Claisen adapters allow an additional entry into single or multiple neck flasks and are useful for reflux with addition reactions.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
ML-70-600	ST 14/10	1	\$57.56



Three-Way Claisen Adapters

These borosilicate glass three-way adapters have two parallel, upper outer Standard Taper joints and one lower inner Standard Taper joint - all of the same size. Claisen adapters allow an additional entry into single or multiple neck flasks and are useful for reflux with addition reactions.

CATALOG NO.	JOINT	OFFSET CENTER TO CENTER	QUANTITY / PACK	PRICE
LG-1550-200	ST 24/40	65mm	1	\$75.96
LG-1550-202	ST 29/42	65mm	1	\$191.34
LG-1550-100	ST 24/40	100mm	1	\$81.29
LG-1550-102	ST 29/42	100mm	1	\$93.92
LG-1550-104	ST 29/26	100mm	1	\$198.97



Full Length Three-Way Thermometer Adapter with 75° Angle Side Arm

With 7mm threaded fitting at top for inert and greaseless thermometer connection, this adapter allows for easy adjustment of thermometer position. Both side and bottom Standard Taper joints are the same size with side joint at 75° angle. Supplied with cap and O-ring.

CATALOG NO.	ACCOMODATION RANGE	APPROX OVERALL HEIGHT	APPROX WIDTH	QUANTITY / PACK	PRICE
LG-1521-100	5/8mm	155mm	100mm	1	\$80.05



90 Degree Bend Adapters with Hose Connection and Standard Taper Joint

These borosilicate glass adapters have a 90° bend to hose connector at top and a Standard Taper inner at bottom.

CATALOG NO.	HOSE CONNECTION O.D.	JOINT	QUANTITY / PACK	PRICE
LG-1702-102	3/8"	ST 24/40	1	\$27.30

POURING ADAPTERS



Pour Out Adapter

Designed for easy transfer of liquids, this borosilicate glass pouring adapter has a Standard Taper 24/40 inner joint to avoid liquid contact with ground zone area of vessels equipped with Standard Taper outer joints.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1353-100	ST 24/40	1	\$35.04



Pouring Adapter

Designed for easy transfer of liquids, this borosilicate glass pouring adapter has a Standard Taper 24/40 inner joint to avoid liquid contact with ground zone area of vessels equipped with Standard Taper outer joints.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1352-100	ST 24/40	1	\$45.06

THREADED ADAPTERS



Threaded Adapters

Featuring a Standard Taper inner joint at bottom and a threaded top, these borosilicate glass adapters connect to threaded nylon bushings which make an O-ring compression seal with thermometers, bleed tubes, etc. when screwed into the adapter. Adapters accommodate thermometers of 7mm O.D. and can be used for vacuum work (for micro sizes see ML-282). Complete set includes glass adapter, nylon bushing (LG-1991-100)*, and Viton® O-ring (LG-1022-502). *Also available with solid bushing (LG-1991-102).

CATALOG NO.	DESCRIPTION	JOINT	QUANTITY / PACK	PRICE
ML-282-800	Complete Adapter Set	ST 14/20	1	\$59.42
ML-282-802	Complete Adapter Set	ST 19/22	1	\$53.53
LG-1986-200	Complete Adapter Set	ST 24/40	1	\$36.71
LG-1986-202	Complete Adapter Set	ST 29/42	1	\$39.96
ML-282-700	Threaded Adapter	ST 14/20	1	\$37.20
LG-1986-100	Threaded Adapter	ST 24/40	1	\$45.58
LG-1986-102	Threaded Adapter	ST 29/42	1	\$48.85
LG-1991-100	Threaded Bushing with Hole		1	\$8.86
LG-1991-102	Threaded Bushing, Solid		1	\$10.15