

LOW FORM BEAKERS



Low Form Graduated Griffin Beakers

Graduated Griffin low form beakers are made of borosilicate glass and have a pour spout.

NOTE: 10L beaker is NOT graduated.

CATALOG NO.	DESCRIPTION	CAPACITY	O.D.	HEIGHT	QUANTITY / PACK	PRICE
LG-3000-100	Graduated Beaker	5mL	22mm	30mm	1	\$8.18
LG-3000-102	Graduated Beaker	5L	170mm	270mm	1	\$94.06
LG-3000-104	10L Beaker - not graduated	10L	217mm	350mm	1	\$252.29



Heavy Duty Low Form Beakers

Low form, heavy duty beakers are made of borosilicate glass and have pour spout.

PLEASE NOTE: 5000mL is graduated; other sizes are not.

CATALOG NO.	CAPACITY	OUTER DIAMETER	HEIGHT	QUANTITY / PACK	PRICE
LG-3002-100	5L	182mm	256mm	1	\$221.15
LG-3002-102	10L	225mm	340mm	1	\$442.33
LG-3002-104	15L	260mm	390mm	1	\$625.58
LG-3002-106	20L	285mm	430mm	1	\$695.08

PTFE BEAKERS



PTFE Beakers

These beakers are molded from PTFE and have a smooth internal finish. The base has been machined to assure a flat contact surface and good heat transfer. With a Griffin profile, they are recommended for trace metal analysis.

CATALOG NO.	CAPACITY	HEIGHT x DIAMETER	QUANTITY / PACK	PRICE
LG-3305-100	25mL	47 x 32mm	1	\$20.25
LG-3305-102	50mL	60 x 43mm	1	\$21.95
LG-3305-104	100mL	68 x 54mm	1	\$28.48
LG-3305-106	150mL	69 x 59mm	1	\$33.78
LG-3305-108	250mL	97 x 66mm	1	\$41.90
LG-3305-110	400mL	106 x 80mm	1	\$69.37
LG-3305-112	500mL	125 x 80mm	1	\$70.68
LG-3305-114	600mL	125 x 90mm	1	\$112.90
LG-3305-116	1000mL	155 x 100mm	1	\$173.75
LG-3305-118	2000mL	210 x 125mm	1	\$272.94
LG-3305-120	5000mL	280 x 170mm	1	\$727.11

JACKETED BEAKERS



Jacketed Beakers

Made of borosilicate glass, these laboratory beakers are fully jacketed to allow temperature control over a wide range. They are similar to a standard Griffin Beaker, but do not have a pouring spout.

CATALOG NO.	CAPACITY	APPROX HEIGHT	APPROX I.D.	HOSE CONNECTION	QUANTITY / PACK	PRICE
LG-3300-100	250mL	95mm	70mm	10mm	1	\$125.33
LG-3300-102	400mL	120mm	75mm	10mm	1	\$144.19
LG-3300-104	1000mL	160mm	110mm	1/2"	1	\$179.10
LG-3300-106	2000mL	195mm	130mm	1/2"	1	\$214.96
LG-3300-108	100mL	75mm	51mm	10mm	1	\$238.99
LG-3300-110	600mL	125mm	90mm	10mm	1	\$406.97