

CRUCIBLE



High Form Gooch Crucibles

With fused-in fritted glass disc, these high form Gooch crucibles, made of borosilicate glass, are for analytical work involving precipitation, drying, and weighing. They are particularly suitable when precipitates will not be dried at temperatures above 110°C and may be used at higher temperatures. Fine (F), Medium (M), and Coarse (C) porosities available. Approximate height 60mm.

Note: Avoid sudden temperature changes when heating or cooling.

CATALOG NO.	CAPACITY	APPROX DIAMETER OF DISC	APPROX HEIGHT	POROSITY	QUANTITY / PACK	PRICE
LG-5470-100	30mL	30mm	60mm	C	1	\$40.21
LG-5470-102	30mL	30mm	60mm	M	1	\$40.21
LG-5470-104	30mL	30mm	60mm	F	1	\$43.81
LG-5470-106	50mL	40mm	60mm	C	1	\$51.68
LG-5470-108	50mL	40mm	60mm	M	1	\$52.89



Low Form Gooch Crucibles

With sealed-in fritted glass disc, these low form Gooch crucibles, made of borosilicate glass, are for analytical work involving precipitation, drying, and weighing. Fine (F), medium (M), and coarse (C) porosities available.

Note: Avoid sudden temperature changes when heating or cooling.

CATALOG NO.	CAPACITY	APPROX DIAMETER OF DISC	APPROX HEIGHT	POROSITY	QUANTITY / PACK	PRICE
LG-5480-102	15mL	20mm	44mm	M	1	\$35.84
LG-5480-106	30mL	30mm	50mm	C	1	\$44.84
LG-5480-108	30mL	30mm	50mm	M	1	\$42.70
LG-5480-110	30mL	30mm	50mm	F	1	\$41.28

FILTER FLASK



Heavy Wall Filter Flasks

Erlenmeyer flask with 3/8" side hose connection and Standard Taper outer joint.

CATALOG NO.	CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-7800-100	125mL	ST 24/40	1	\$85.43
LG-7800-102	250mL	ST 24/40	1	\$95.45
LG-7800-104	500mL	ST 24/40	1	\$109.76
LG-7800-106	1000mL	ST 24/40	1	\$142.38
LG-7800-110	500mL	ST 29/42	1	\$107.46
LG-7800-112	1000mL	ST 29/42	1	\$148.70
LG-7800-114	2000mL	ST 29/42	1	\$192.93



Duran Filter Flasks with Hose Connection

Duran heavy wall filter flask with removable hose connection. Also available safety coated.

CATALOG NO.	DESCRIPTION	CAPACITY	I.D.	O.D.	HEIGHT	QUANTITY / PACK	PRICE
LG-7802-096	Non-coated Flask	3000mL	60mm	170mm	295mm	1	\$189.27
LG-7802-098	Non-coated Flask	5000mL	70mm	185mm	360mm	1	\$293.91
LG-7802-100	Non-coated Flask	10000mL	70mm	240mm	420mm	1	\$447.80
LG-7802-102	Non-coated Flask	15000mL	70mm	255mm	500mm	1	\$479.77
LG-7802-104	Non-coated Flask	20000mL	70mm	290mm	535mm	1	\$583.31
LG-7802-150	Replacement Hose Connection Set					1	\$6.84
LG-7802-196	Coated Flask	3000mL	60mm	170mm	295mm	1	\$240.16
LG-7802-198	Coated Flask	5000mL	70mm	185mm	360mm	1	\$350.17
LG-7802-200	Coated Flask	10000mL	70mm	240mm	420mm	1	\$498.29
LG-7802-202	Coated Flask	15000mL	70mm	255mm	500mm	1	\$567.74
LG-7802-204	Coated Flask	20000mL	70mm	290mm	535mm	1	\$675.21

FILTERWARE Filter Flasks with Stopper Joint

These borosilicate glass heavy wall filtering flasks can be used with FILTERWARE® Filtration Apparatus (LG-1755). Supplied with a silicone stopper with 9/16" hole and a #2 hose connection, the 1000mL flask has a 3/8" hose connection, whereas the 2000mL and 4000mL flasks have 1/2" connections.



CATALOG NO.	CAPACITY	STOPPER SIZE	QUANTITY / PACK	PRICE
BP-1760-000	1000mL	#8 - Silicone	1	\$73.99
BP-1760-002	2000mL	#9 - Silicone	1	\$126.98
BP-1760-004	4000mL	#12 - Silicone	1	\$183.60

FILTER FUNNEL




Buchner Filter Funnels with Fritted Disc

Used in applications where filter paper would be attacked by the material being filtered, these borosilicate glass Buchner filter funnels with sealed-in fritted glass discs are available in various capacities* with Coarse(C), Medium(M), or Fine(F) porosities.

*Larger sizes available upon request.

CATALOG NO.	CAPACITY	APPROX DISC DIAMETER	APPROX HEIGHT ABOVE DISC	POROSITY	QUANTITY / PACK	PRICE
LG-7080-102	2mL	10mm	30mm	Coarse	1	\$44.97
LG-7080-104	15mL	20mm	45mm	Coarse	1	\$52.44
LG-7080-106	30mL	30mm	45mm	Coarse	1	\$59.17
LG-7080-108	60mL	40mm	50mm	Coarse	1	\$62.66
LG-7080-110	150mL	60mm	60mm	Coarse	1	\$114.38
LG-7080-112	350mL	80mm	80mm	Coarse	1	\$154.39
LG-7080-114	600mL	90mm	90mm	Coarse	1	\$179.51
LG-7080-116	2000mL	125mm	147mm	Coarse	1	\$409.39
LG-7080-118	3000mL	150mm	157mm	Coarse	1	\$568.10
LG-7080-132	2mL	10mm	30mm	Medium	1	\$44.97
LG-7080-134	15mL	20mm	45mm	Medium	1	\$51.24
LG-7080-136	30mL	30mm	45mm	Medium	1	\$57.82
LG-7080-138	60mL	40mm	50mm	Medium	1	\$61.54
LG-7080-140	150mL	60mm	60mm	Medium	1	\$112.77
LG-7080-142	350mL	80mm	80mm	Medium	1	\$154.39
LG-7080-144	600mL	90mm	90mm	Medium	1	\$179.51
LG-7080-146	2000mL	125mm	147mm	Medium	1	\$409.39
LG-7080-148	3000mL	150mm	157mm	Medium	1	\$568.10
LG-7080-162	2mL	10mm	30mm	Fine	1	\$44.97
LG-7080-164	15mL	20mm	45mm	Fine	1	\$51.24
LG-7080-166	30mL	30mm	45mm	Fine	1	\$59.17
LG-7080-168	60mL	40mm	50mm	Fine	1	\$63.77
LG-7080-170	150mL	60mm	60mm	Fine	1	\$116.54
LG-7080-172	350mL	80mm	80mm	Fine	1	\$146.90
LG-7080-174	600mL	90mm	90mm	Fine	1	\$182.41
LG-7080-176	2000mL	125mm	147mm	Fine	1	\$497.77
LG-7080-178	3000mL	150mm	157mm	Fine	1	\$596.36



Jacketed Buchner Filter Funnels with Fritted Disc

The jacket of these Buchner filter funnels encloses entire funnel bowl and upper part of stem. Circulation of liquid through the jacket permits warming or cooling of funnel and contents during filtration. Made of borosilicate glass with sealed-in fritted glass disc, they are available in coarse (C), medium (M) and fine (F) porosities.

CATALOG NO.	CAPACITY	APPROX DIAMETER OF DISC	APPROX HEIGHT ABOVE DISC	HOSE CONNECTION	POROSITY	QUANTITY / PACK	PRICE
LG-7110-100	30mL	30mm	45mm	3/8"	C	1	\$342.78
LG-7110-102	60mL	40mm	50mm	3/8"	C	1	\$399.02
LG-7110-106	600mL	90mm	90mm	1/2"	C	1	\$568.65
LG-7110-108	2000mL	125mm	147mm	1/2"	C	1	\$988.72
LG-7110-112	60mL	40mm	50mm	3/8"	M	1	\$399.02
LG-7110-114	150mL	60mm	50mm	3/8"	M	1	\$425.11
LG-7110-116	600mL	90mm	90mm	1/2"	M	1	\$568.64
LG-7110-126	600mL	90mm	90mm	1/2"	F	1	\$630.25
LG-7110-128	2000mL	125mm	147mm	1/2"	F	1	\$988.72



Buchner Filter Funnels with Fritted Disc and Vacuum Adapter

These Buchner filter funnels with vacuum adapter feature a full length Standard Taper inner joint sealed to the stem with side hose connection just above the grinding. Made of borosilicate glass with sealed-in fritted glass disc, they are available in coarse (C) and medium (M) porosities. Recommended for use with LG-7820 heavy wall filter flasks.

CATALOG NO.	CAPACITY	APPROX DIAMETER OF DISC	POROSITY	JOINT	QUANTITY / PACK	PRICE
LG-7120-098	15mL	20mm	C	ST 29/26	1	\$148.67
LG-7120-100	15mL	20mm	C	ST 24/40	1	\$158.78
LG-7120-101	30mL	30mm	C	ST 29/26	1	\$164.95
LG-7120-102	30mL	30mm	C	ST 24/40	1	\$168.51
LG-7120-103	60mL	40mm	C	ST 24/40	1	\$164.74
LG-7120-104	60mL	40mm	C	ST 29/42	1	\$134.89
LG-7120-105	150mL	60mm	C	ST 24/40	1	\$217.23
LG-7120-106	150mL	60mm	C	ST 29/42	1	\$214.55
LG-7120-107	350mL	80mm	C	ST 24/40	1	\$270.68
LG-7120-108	350mL	80mm	C	ST 29/42	1	\$261.76
LG-7120-109	600mL	90mm	C	ST 24/40	1	\$305.09
LG-7120-110	600mL	90mm	C	ST 29/42	1	\$314.07
LG-7120-111	60mL	40mm	C	ST 29/26	1	\$174.31
LG-7120-112	150mL	60mm	C	ST 29/26	1	\$224.19
LG-7120-118	15mL	20mm	M	ST 29/26	1	\$148.67
LG-7120-120	15mL	20mm	M	ST 24/40	1	\$158.78
LG-7120-121	30mL	30mm	M	ST 29/26	1	\$164.95
LG-7120-122	30mL	30mm	M	ST 24/40	1	\$168.51
LG-7120-123	60mL	40mm	M	ST 24/40	1	\$170.37
LG-7120-124	60mL	40mm	M	ST 29/42	1	\$166.95
LG-7120-125	150mL	60mm	M	ST 24/40	1	\$227.69
LG-7120-126	150mL	60mm	M	ST 29/42	1	\$224.19
LG-7120-127	350mL	80mm	M	ST 24/40	1	\$270.68
LG-7120-129	600mL	90mm	M	ST 24/40	1	\$305.09
LG-7120-130	600mL	90mm	M	ST 29/42	1	\$314.09
LG-7120-131	60mL	40mm	M	ST 29/26	1	\$174.31
LG-7120-132	150mL	60mm	M	ST 29/26	1	\$224.19
LG-7120-133	350mL	80mm	M	ST 29/26	1	\$261.76
LG-7120-134	600mL	90mm	M	ST 29/26	1	\$314.09
ML-1060-700	15mL	20mm	C	ST 14/20	1	\$102.38
ML-1060-704	30mL	30mm	C	ST 14/20	1	\$107.74
ML-1060-708	60mL	40mm	C	ST 14/20	1	\$115.67
ML-1060-709	60mL	40mm	C	ST 19/22	1	\$119.61
ML-1060-710	15mL	20mm	M	ST 14/20	1	\$102.38
ML-1060-712	15mL	20mm	M	ST 19/22	1	\$106.93
ML-1060-714	30mL	30mm	M	ST 14/20	1	\$112.98
ML-1060-718	60mL	40mm	M	ST 14/20	1	\$119.17

FILTER TUBE



Straight Filter Tubes

Made of borosilicate glass, these filtration tubes are straight, 200mm long, and have disc sealed-in at center. They are intended for sealing to apparatus requiring filtration. Tube should not be heated too close to disc.

CATALOG NO.	O.D.	POROSITY	APPROX DIAMETER OF DISC	QUANTITY / PACK	PRICE
LG-7160-120	13mm	Coarse	10mm	1	\$33.70
LG-7160-122	25mm	Coarse	20mm	1	\$39.47
LG-7160-124	35mm	Coarse	30mm	1	\$44.55
LG-7160-140	13mm	Medium	10mm	1	\$33.70
LG-7160-142	25mm	Medium	20mm	1	\$39.47
LG-7160-160	13mm	Fine	10mm	1	\$33.70
LG-7160-162	25mm	Fine	20mm	1	\$39.47
LG-7160-164	35mm	Fine	30mm	1	\$44.55



Filter Tubes with Fritted Disc

These gas dispersion filter tubes consist of a closed bottom fritted glass cylinder. Made of borosilicate glass, they are 5mm in diameter and 135mm long. Useful when gas must be introduced into a container with narrow opening.

CATALOG NO.	POROSITY	QUANTITY / PACK	PRICE
ML-1080-700	Fine	1	\$27.98
ML-1080-702	Medium	1	\$30.76
ML-1080-704	Coarse	1	\$27.98



Filter Tubes with 10mm Fritted Disc

These filter tubes, made of borosilicate glass, consist of a fritted disc sealed at the end of a 6mm O.D. tube. They are useful in removing supernatant liquid from above a precipitate. Available in coarse (C), medium (M) and fine (F) porosities.

CATALOG NO.	POROSITY	QUANTITY / PACK	PRICE
ML-1090-700	Coarse	1	\$49.30
ML-1090-710	Medium	1	\$44.33
ML-1090-720	Fine	1	\$50.23

FILTERWARE



FILTERWARE Funnel Apparatus

These filter assemblies for FILTERWARE™ filtration system consist of an all-borosilicate glass funnel, fritted glass support base, aluminum clamp, and silicone stopper. Choose from 47mm or 90mm filter assembly sizes. Individual components are offered in various sizes for a customized unit. For use with BP-1827 filter flasks, sold separately.

NOTE: Aluminum Clamps BP-1753-000, BP-1753-090 and Graduated Funnels BP-1751-000, BP-1751-090 are also components of BP-1825 FILTERWARE™ All-Glass Filtration Apparatus.

CATALOG NO.	DESCRIPTION	CAPACITY	LENGTH	JOINT SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
BP-1755-000	Complete Filter Apparatus	300mL	47mm	40/35		1	\$252.11
BP-1755-090	Complete Filter Apparatus	1000mL	90mm	40/35		1	\$649.45
BP-1751-000	Graduated Funnel	300mL	47mm			1	\$78.12
BP-1781-090	Graduated Funnel	1000mL	90mm			1	\$172.16
BP-1752-001	Fritted Glass Support Base, 47mm Filter Apparatus			40/35		1	\$118.77
BP-1752-091	Fritted Glass Support Base, 90mm Filter Apparatus			40/35		1	\$336.54
BP-1753-000	Aluminum Clamp, 47mm Filter Apparatus					1	\$78.80
BP-1753-090	Aluminum Clamp, 90mm Filter Apparatus					1	\$111.74
BP-1763-000	Silicone Stopper for 1000mL Flask, 9/16" hole				8	1	\$21.71
BP-1763-002	Silicone Stopper for 2000mL Flask, 9/16" hole				9	1	\$22.97
BP-1763-004	Silicone Stopper for 4000mL Flask, 9/16" hole				12	1	\$27.53



FILTERWARE All-Glass Filter Apparatus

Choose from 47mm or 90mm microfiltration FILTERWARE™ all-borosilicate glass filter apparatus consisting of a graduated funnel, a Standard Taper 40/35 glass joint connecting the coarse fritted filter support base to the filter flask, and an aluminum clamp (BP-1753). Use for HPLC mobile phase filtration, analysis of particulate or bacterial contamination, or for general laboratory microfiltration.

NOTE: Components of 47mm apparatus are interchangeable with the corresponding components of other manufacturer's items and can be used with filtering manifolds.

CATALOG NO.	DESCRIPTION	FLASK CAPACITY	FUNNEL CAPACITY	JOINT	QUANTITY / PACK	PRICE
BP-1825-000	47mm Filter Apparatus	1L	300mL		1	\$445.10
BP-1825-002	47mm Filter Apparatus	2L	500mL		1	\$577.55
BP-1825-004	47mm Filter Apparatus	4L	1000mL		1	\$630.55
BP-1825-090	90mm Filter Apparatus	4L	1000mL		1	\$949.31
BP-1771-000	Filter Funnel, 47mm Apparatus		500mL		1	\$122.33
BP-1781-000	Filter Funnel, 47mm Apparatus		1000mL		1	\$141.46
BP-1781-092	Filter Funnel, 90mm Apparatus		2000mL		1	\$315.74
BP-1826-000	Coarse Fritted Glass Support Base, 47mm Apparatus			ST 40/35	1	\$205.64
BP-1826-090	Coarse Fritted Glass Support Base, 90mm Apparatus			ST 40/35	1	\$547.76
BP-1827-000	Filter Flask	1000mL		ST 40/35	1	\$140.76
BP-1827-002	Filter Flask	2000mL		ST 40/35	1	\$218.77
BP-1827-004	Filter Flask	4000mL		ST 40/35	1	\$249.73
BP-1828-000	Filter Funnel, 47mm Apparatus		2000mL		1	\$229.62

VACUUM TAKE-OFF



Hollow Glass Stoppers

CATALOG NO.	STOPPER NUMBER	QUANTITY / PACK	PRICE
LG-10326-108	22	1	\$22.88
LG-10326-110	27	1	\$27.45

VACUUM TAKE-OFF



Straight Vacuum Adapters

For vacuum filtration into vessels. 3/8" hose connector is provided for attachment to vacuum source. Top opening accommodates a #6 stopper.

CATALOG NO.	JOINT	STOPPER NUMBER	QUANTITY / PACK	PRICE
LG-1272-100	ST 24/40	6	1	\$70.96
LG-1272-102	ST 29/42	6	1	\$74.81



Filtration Vacuum Adapters

Adapter for reduced pressure filtration using plain stem Buchner funnels. Top is tooled to accept special neoprene collar for funnels up to 4000mL capacity. Funnel sets in collar for vacuum seal and is easy to remove. See LG-1275C for collars.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1275-100	ST 24/40	3/8"	1	\$88.27
LG-1275-102	ST 29/42	3/8"	1	\$93.03
ML-265-700	ST 14/20	5/16"	1	\$82.62
ML-265-702	ST 19/22	5/16"	1	\$85.85



Straight Vacuum Adapters with 3/8 Inch Hose Connection

With a standard taper outer joint at one end and standard taper inner joint at the other; hose connection for vacuum source.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1330-100	ST 24/40 Inner, ST 24/40 Outer	3/8"	1	\$60.87
LG-1330-102	ST 29/42 Inner, ST 29/42 Outer	3/8"	1	\$67.09


105 Degree Bent Vacuum Adapter



Adapter with a drip tip and a 3/8" hose connection.

CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
LG-1371-100	ST 24/40	1	\$103.73

Straight Vacuum Adapters with Hose Connection and Sealed Delivery Tube



With a Standard Taper outer joint at top and a Standard Taper inner joint of the same size at bottom, 3/8" side hose connection and sealed delivery tube.


CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1840-100	ST 24/40 Outer, ST 24/40 Inner	3/8"	1	\$74.35
LG-1840-102	ST 29/42 Outer, ST 29/42 Inner	3/8"	1	\$80.43
ML-375-700	ST 14/20 Outer, ST 14/20 Inner	3/8"	1	\$58.28
ML-375-702	ST 19/22 Outer, ST 19/22 Inner	3/8"	1	\$54.10

Straight Vacuum Adapters with 5/16 Inch Hose Connection



Sealed in delivery tube extends 90mm below lower inner joint. Upper and lower joints are the same. Supplied with 5/16" hose connection.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-1842-100	ST 24/40	5/16"	1	\$82.04
LG-1842-102	ST 29/42	5/16"	1	\$78.83
ML-380-700	ST 14/20	5/16"	1	\$58.04
ML-380-702	ST 19/22	5/16"	1	\$63.46



90 Degree Bent Gas Inlet Vacuum Adapters

Provided with a standard taper inner joint and delivery tube extending 250mm below lower joint. 3/8" hose connectors allow attachment of flexible tubing.

CATALOG NO.	JOINT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-2050-102	ST 24/40	3/8"	1	\$108.34
LG-2050-104	ST 29/42	3/8"	1	\$116.07

STOPCOCK



Balloon Adapters

These Air-Tite borosilicate glass adapters are designed with a T-bore stopcock for addition or venting inert atmospheres by attaching a balloon to one of two stems. Choose from different sized Standard Taper joints and plug type. IG-22 has a 4mm glass stopcock while IG-23T has a 4mm PTFE plug.

CATALOG NO.	DESCRIPTION	JOINT	QUANTITY / PACK	PRICE
IG-22-100	Glass Stopcock	ST 14/20	1	\$147.79
IG-22-104	Glass Stopcock	ST 24/40	1	\$159.45
IG-23T-100	PTFE Plug	ST 14/20	1	\$151.03
IG-23T-104	PTFE Plug	ST 24/40	1	\$157.89



Adapters with 90 Degree Bend to Hose Connection and Stopcock

Made of borosilicate glass, these adapters have a 90° bend to 5/16" O.D. hose connection on top and inner Standard Taper joint on bottom. Equipped with a stopcock, they are used to control gas/liquid addition in reaction systems. ML-210 and LG-1710 have glass plugs; ML-211T and LG-1710T have PTFE plugs.

CATALOG NO.	JOINT	PLUG TYPE	QUANTITY / PACK	PRICE
LG-1710-102	ST 24/40	Glass	1	\$73.04
LG-1710-104	ST 29/42	Glass	1	\$82.62
LG-1710T-202	ST 24/40	PTFE	1	\$99.46
LG-1710T-204	ST 29/42	PTFE	1	\$103.34
ML-210-700	ST 14/20	Glass	1	\$68.80
ML-210-702	ST 19/22	Glass	1	\$66.79
ML-211T-700	ST 14/20	PTFE	1	\$89.98
ML-211T-702	ST 19/22	PTFE	1	\$98.95

Air-Tite Adapters with Vacuum Valve

These borosilicate glass adapters are equipped with vacuum valve for closing joints yet allowing for reconnection to other equipment with minimum exposure to air and vapors. The PTFE valve minimizes contamination and maintains a vacuum in the micron region. Standard Taper upper outer and bottom inner are same size.



CATALOG NO.	JOINT	QUANTITY / PACK	PRICE
IG-202-100	ST 14/20	1	\$183.21
IG-202-102	ST 19/22	1	\$174.09
IG-202-104	ST 24/40	1	\$188.66

FRITTED GLASS DISK

Fritted Discs

Used in manufacturing of laboratory glassware, these borosilicate glass fritted discs are available in four porosities - Fine (4/5.5µm), Medium (10/15µm), Coarse (40/60µm), and Extra Coarse (170/220µm).



CATALOG NO.	APPROX O.D.	POROSITY	QUANTITY / PACK	PRICE
LG-8675-118	40mm	Extra Coarse	1	\$13.77
LG-8675-122	60mm	Extra Coarse	1	\$23.09
LG-8675-140	5mm	Coarse	1	\$10.56
LG-8675-152	20mm	Coarse	1	\$9.52
LG-8675-156	30mm	Coarse	1	\$16.30
LG-8675-158	40mm	Coarse	1	\$24.83
LG-8675-162	60mm	Coarse	1	\$23.62
LG-8675-180	5mm	Medium	1	\$9.44
LG-8675-192	20mm	Medium	1	\$9.62
LG-8675-196	30mm	Medium	1	\$10.35
LG-8675-198	40mm	Medium	1	\$27.32
LG-8675-202	60mm	Medium	1	\$23.64


O-RING CONNECTOR

Air-Tite Recovery Flasks with O-Ring Connector

Made of borosilicate glass, this recovery style flask has a size 15 O-ring joint sealed on the opening. Supplied with LG-1022-530 Viton® O-ring.




CATALOG NO.	CAPACITY	QUANTITY / PACK	PRICE
LG-11081-100	10mL	1	\$34.08
LG-11081-102	25mL	1	\$34.08
LG-11081-104	50mL	1	\$34.08
LG-11081-106	100mL	1	\$36.52



Air-Tite Connecting/Inlet Adapter with Vacuum Valve

This borosilicate glass adapter has a right angle configuration with a size 15 O-ring joint sealed below valve stem and a side inlet. It is easily connected to vacuum manifolds and supplied with LG-1022-530 Viton® O-ring.

CATALOG NO.	APPROXIMATE OVERALL DIMENSIONS (H x W)	APPROXIMATE SIDEARM O.D.	VALVE SIZE	QUANTITY / PACK	PRICE
LG-11078-100	100 x 58mm	9mm	0-4mm	1	\$108.32



Air-Tite Two O-Ring Connecting Adapter with Vacuum Valve

With a right angle configuration, this borosilicate glass connecting adapter features two size 15 O-ring joints below valve stem for completing a seal to vacuum manifold or related apparatus. Supplied with two LG-1022-530 Viton O-rings.

CATALOG NO.	APPROXIMATE OVERALL DIMENSIONS (H x W)	VALVE SIZE	QUANTITY / PACK	PRICE
LG-11079-100	100 x 60mm	0-4mm	1	\$100.83

O-RING CONNECTOR



Air-Tite Hose Connector Adapter

This borosilicate glass adapter is supplied as a connection to LG-11071 and LG-11072 manifolds. Offered with size 15 O-ring joint for connection to heavy wall rubber tubing used with vacuum pumps and traps. Supplied with a LG-1022-530 Viton O-ring.

CATALOG NO.	APPROX OVERALL HEIGHT	HOSE CONNECTION O.D.	QUANTITY / PACK	PRICE
LG-11073-100	50mm	3/8"	1	\$23.65

O-RING CONNECTOR



Air-Tite Connecting Adapter with O-Ring

Made of borosilicate glass, this connecting adapter has a Standard Taper inner joint on one end for connection to apparatus with outer joint and a size 15 O-ring joint on other end for connection to vacuum manifolds. Supplied with one LG-1022-530 Viton O-ring.

CATALOG NO.	APPROXIMATE OVERALL HEIGHT	JOINT	QUANTITY / PACK	PRICE
LG-11074-100	60mm	ST 14/20	1	\$38.59



Air-Tite Connecting Adapters with O-Ring

Made of borosilicate glass, these connecting adapters have Standard Taper outer joint on one end for connection to apparatus with inner joint and a size 15 O-ring joint on other end for connection to vacuum manifolds. Supplied with one LG-1022-530 Viton O-ring.

CATALOG NO.	APPROXIMATE OVERALL HEIGHT	JOINT	QUANTITY / PACK	PRICE
LG-11075-100	60mm	ST 14/20	1	\$39.65
LG-11075-102	90mm	ST 24/40	1	\$40.88

O-RING CONNECTOR



Glass Stopper with Size 15 O-Ring

Size 15 O-ring joint is closed and provides a seal for open ports on the manifold. Stopper can be used as a closure on the in-line filter or similar opening flasks. Supplied with one LG-1022-530 Viton® O-ring.

CATALOG NO.	OVERALL HEIGHT	O-RING SIZE	QUANTITY / PACK	PRICE
LG-11082-100	25mm	15	1	\$23.49