

ADDITION FUNNELS



Cylindrical Addition Funnels with Pressure Equalizing Sidearm

Featuring a PTFE stopcock, these borosilicate glass cylindrical funnels have a pressure equalizing sidearm. Standard Taper outer joint at the top is the same size as the lower Standard Taper inner joint. Inner joint has a sealed-in drip tube.

CATALOG NO.	CAPACITY	STOPCOCK SIZE	JOINT SIZE	QUANTITY / PACK	PRICE
LG-8266T-100	60mL	PTFE 2	ST 24/40	1	\$252.59
LG-8266T-102	125mL	PTFE 2	ST 24/40	1	\$255.08
LG-8266T-104	250mL	PTFE 2	ST 24/40	1	\$275.15
LG-8266T-106	500mL	PTFE 4	ST 24/40	1	\$298.91
LG-8266T-107	500mL	PTFE 4	ST 29/26	1	\$310.43
LG-8266T-108	500mL	PTFE 4	ST 29/42	1	\$310.38
LG-8266T-110	1000mL	PTFE 4	ST 24/40	1	\$348.25
LG-8266T-111	1000mL	PTFE 4	ST 29/26	1	\$359.37
LG-8266T-112	1000mL	PTFE 4	ST 29/42	1	\$359.32
ML-1421T-700	10mL	PTFE 2	ST 14/20	1	\$189.79
ML-1421T-702	25mL	PTFE 2	ST 14/20	1	\$189.79
ML-1421T-704	50mL	PTFE 2	ST 14/20	1	\$189.79
ML-1421T-706	60mL	PTFE 2	ST 14/20	1	\$189.79
ML-1421T-750	10mL	PTFE 2	ST 19/22	1	\$183.08
ML-1421T-752	25mL	PTFE 2	ST 19/22	1	\$183.08
ML-1421T-756	60mL	PTFE 2	ST 19/22	1	\$183.08



Graduated Cylindrical Addition Funnel with Pressure Equalizing Sidearm

Featuring graduations with a double scale from the PTFE stopcock, these borosilicate glass cylindrical funnels have a pressure equalizing sidearm. Standard Taper outer joint at the top is the same size as the lower Standard Taper inner joint.

CATALOG NO.	CAPACITY	GRADUATION	STOPCOCK SIZE	JOINT SIZE	QUANTITY / PACK	PRICE
ML-1431T-700	10mL	0.2mL	PTFE 2	ST 14/20	1	\$231.42
ML-1431T-750	10mL	0.2mL	PTFE 2	ST 19/22	1	\$238.74
ML-1431T-702	25mL	0.5mL	PTFE 2	ST 14/20	1	\$231.40
ML-1431T-752	25mL	0.5mL	PTFE 2	ST 19/22	1	\$254.97
ML-1431T-754	50mL	1mL	PTFE 2	ST 19/22	1	\$250.92
ML-1431T-704	50mL	1.0mL	PTFE 2	ST 14/20	1	\$246.49
ML-1431T-756	60mL	1mL	PTFE 2	ST 19/22	1	\$254.97
ML-1431T-706	60mL	1mL	PTFE 2	ST 14/20	1	\$250.48
LG-8271T-100	60mL	1.0mL	PTFE 2	ST 24/40	1	\$264.41
LG-8271T-102	125mL	1.0mL	PTFE 2	ST 24/40	1	\$268.37
LG-8271T-104	250mL	5.0mL	PTFE 2	ST 24/40	1	\$305.94
LG-8271T-106	500mL	5.0mL	PTFE 4	ST 24/40	1	\$364.03
LG-8271T-108	500mL	5.0mL	PTFE 4	ST 29/42	1	\$375.52
LG-8271T-110	1000mL	10mL	PTFE 4	ST 24/40	1	\$437.08



Pressure Equalizing Cylindrical Addition Funnel with Two Stopcocks, PTFE Plugs

These cylindrical addition funnels have a pressure equalizing sidearm with a Standard Taper 2mm PTFE stopcock, allowing the funnel to be refilled without opening the system to the atmosphere. They feature the same size Standard Taper joint at top and bottom and a lower stopcock with PTFE plug and sealed-in drip tube.

CATALOG NO.	CAPACITY	JOINT	LOWER STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8268T-100	60mL	ST 24/40	PTFE 2mm	1	\$310.19
LG-8268T-102	125mL	ST 24/40	PTFE 2mm	1	\$314.59
LG-8268T-104	250mL	ST 24/40	PTFE 2mm	1	\$332.53
LG-8268T-106	500mL	ST 24/40	PTFE 4mm	1	\$359.32



Graduated Pressure Equalizing Cylindrical Addition Funnel with Two PTFE Plugs

Graduated with double scale from the bottom PTFE stopcock, these cylindrical addition funnels have a pressure equalizing sidearm with Standard Taper 2mm PTFE stopcock, allowing the funnel to be refilled without opening the system to the atmosphere. They feature the same size Standard Taper joint at top and bottom - bottom with sealed-in drip tube.

CATALOG NO.	CAPACITY	SUBDIVISIONS	JOINT	LOWER STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8273T-100	60mL	1mL	ST 24/40	PTFE 2mm	1	\$327.91
LG-8273T-102	125mL	1mL	ST 24/40	PTFE 2mm	1	\$331.12
LG-8273T-104	250mL	5mL	ST 24/40	PTFE 2mm	1	\$400.34
LG-8273T-106	500mL	5mL	ST 24/40	PTFE 4mm	1	\$432.15
LG-8273T-108	500mL	5mL	ST 29/42	PTFE 4mm	1	\$443.62
LG-8273T-110	1000mL	10mL	ST 24/40	PTFE 4mm	1	\$495.57
LG-8273T-112	1000mL	10mL	ST 29/42	PTFE 4mm	1	\$506.66



Graduated Pressure Equalizing Cylindrical Addition Funnel with PTFE Valve

Graduated with double scale from bottom valve, these cylindrical addition funnels have a pressure equalizing sidearm and a 0-4mm PTFE valve to allow for semi-needle valve action. They feature the same size Standard Taper joint at top and bottom - bottom with sealed-in drip tube.

CATALOG NO.	CAPACITY	SUBDIVISIONS	JOINT	QUANTITY / PACK	PRICE
LG-8269T-100	250mL	5mL	ST 24/40	1	\$313.16
LG-8269T-102	250mL	5mL	ST 29/42	1	\$325.85
LG-8269T-104	500mL	5mL	ST 24/40	1	\$353.06
LG-8269T-106	500mL	5mL	ST 29/42	1	\$365.07
LG-8269T-108	1000mL	10mL	ST 24/40	1	\$414.52
LG-8269T-110	1000mL	10mL	ST 29/42	1	\$422.37

POURING FUNNELS



60° Powder Funnels

These 60° funnels are used for the introduction of liquids or powders into containers with the same Standard Taper joint. The funnel design promotes free flow and provides smooth surface for transfer of material.

CATALOG NO.	APPROX DIAMETER AT TOP	JOINT	QUANTITY / PACK	PRICE
LG-8305-100	75mm	ST 24/40	1	\$50.87
LG-8305-102	100mm	ST 24/40	1	\$55.33
LG-8305-104	100mm	ST 29/42	1	\$63.51
LG-8305-106	125mm	ST 24/40	1	\$75.53
LG-8305-108	125mm	ST 29/42	1	\$74.69
LG-8305-110	150mm	ST 45/50	1	\$93.71
ML-1360-700	50mm	ST 14/20	1	\$42.08



Straight Powder Funnels

These straight funnels are used for the introduction of liquids or powders into containers with the same Standard Taper joint. The funnel design promotes free flow and provides smooth surface for transfer of material.

CATALOG NO.	APPROX DIAMETER AT TOP	JOINT	QUANTITY / PACK	PRICE
LG-8300-102	100mm	ST 24/40	1	\$57.40
LG-8300-103	100mm	ST 29/26	1	\$56.15
LG-8300-105	100mm	ST 45/50	1	\$63.64
LG-8300-106	125mm	ST 24/40	1	\$53.93
LG-8300-108	125mm	ST 29/42	1	\$63.64
LG-8300-110	150mm	ST 45/50	1	\$92.54

PTFE Powder Funnels



These PTFE funnels are totally chemically inert with a super smooth non-stick internal finish. They are suitable for powders and liquids.

CATALOG NO.	DIAMETER	HEIGHT	STEM LENGTH	QUANTITY / PACK	PRICE
LG-8313-100	30mm	50mm	8mm	1	\$19.26
LG-8313-102	50mm	85mm	10mm	1	\$36.80
LG-8313-104	70mm	110mm	12mm	1	\$62.81
LG-8313-106	100mm	170mm	14mm	1	\$114.88
LG-8313-108	150mm	220mm	16mm	1	\$229.77

60° Angled Powder Funnels



These 60° angled funnels are designed for use with multi-neck flasks. Suitable for powders and liquids, they are used to introduce materials into containers with the same Standard Taper joint.

CATALOG NO.	DIAMETER	JOINT	QUANTITY / PACK	PRICE
LG-8312-102	75mm	ST 24/40	1	\$82.63
LG-8312-104	100mm	ST 24/40	1	\$89.59
LG-8312-106	100mm	ST 29/42	1	\$86.77
LG-8312-108	125mm	ST 24/40	1	\$93.47
LG-8312-110	125mm	ST 29/42	1	\$101.73

Offset Powder Funnels



These offset funnels are designed with a straight side and are suggested for use with multi-neck flasks. Suitable for powders and liquids, they are used to introduce materials into containers with the same size Standard Taper joint.

CATALOG NO.	DIAMETER	JOINT	QUANTITY / PACK	PRICE
LG-8303-100	75mm	ST 24/40	1	\$63.86
LG-8303-102	100mm	ST 24/40	1	\$71.77
LG-8303-103	100mm	ST 29/26	1	\$72.35
LG-8303-104	100mm	ST 29/42	1	\$72.35
ML-1365-700	50mm	ST 14/20	1	\$51.47
ML-1365-702	50mm	ST 19/22	1	\$58.86



60° Angled Offset Powder Funnel

These 60° offset funnels are angled for use with multi-neck flasks. They feature a uniform wall and a strong beaded rim. Suitable for powders and liquids, they are used to introduce materials into containers with the same size Standard Taper joint.

CATALOG NO.	DIAMETER	JOINT	QUANTITY / PACK	PRICE
LG-8307-100	75mm	ST 24/40	1	\$66.03
LG-8307-102	100mm	ST 24/40	1	\$79.87
LG-8307-104	100mm	ST 29/42	1	\$92.39

SEPARATORY FUNNELS



Squibb Separatory Funnel with Glass Stopper and PTFE Plug, Drip Tip

Featuring bottom stopcock with PTFE plug and glass stopper, these Squibb separatory funnels have a lower Standard Taper inner joint with a sealed-in drip tube.

CATALOG NO.	CAPACITY	LOWER JOINT	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8431T-102	125mL	ST 24/40	ST 22	PTFE 2mm	1	\$169.20
LG-8431T-104	250mL	ST 24/40	ST 22	PTFE 2mm	1	\$186.29
LG-8431T-106	500mL	ST 24/40	ST 27	PTFE 4mm	1	\$226.71
LG-8431T-108	1000mL	ST 29/42	ST 27	PTFE 4mm	1	\$321.91
LG-8431T-110	2000mL	ST 29/42	ST 38	PTFE 6mm	1	\$423.36
ML-1451T-706	30mL	ST 14/20	ST 14/20	PTFE 2mm	1	\$123.42
ML-1451T-710	60mL	ST 14/20	ST 14/20	PTFE 2mm	1	\$133.28
ML-1451T-762	125mL	ST 19/22	ST 14/20	PTFE 2mm	1	\$128.00



Squibb Separatory Funnels

These economy style Squibb separatory funnels come complete with a detachable PTFE stopcock, glass drip tip, and plastic stopper.

NOTE: Glass tip (LG-3958-106) and PTFE valve (3958T-104) must be purchased together for proper fit.

CATALOG NO.	DESCRIPTION	CAPACITY	STOPPER SIZE	QUANTITY / PACK	PRICE
LG-8435T-100	Complete Separatory Funnel	60mL	ST 16	1	\$95.82
LG-8435T-102	Complete Separatory Funnel	125mL	ST 22	1	\$119.25
LG-8435T-104	Complete Separatory Funnel	250mL	ST 22	1	\$136.27
LG-8435T-106	Complete Separatory Funnel	500mL	ST 27	1	\$147.11
LG-3958T-104	PTFE Stopcock, 2mm			1	\$43.16
LG-8436-100	Replacement Funnel	60mL	ST 16	1	\$53.25
LG-8436-102	Replacement Funnel	125mL	ST 22	1	\$65.99
LG-8436-104	Replacement Funnel	250mL	ST 22	1	\$75.97
LG-8436-106	Replacement Funnel	500mL	ST 27	1	\$87.04



Squibb Separatory Funnels with PTFE Plug and Plastic Stopper

Featuring a bottom stopcock with PTFE plug, these Squibb separatory funnels come with a plastic stopper.

CATALOG NO.	CAPACITY	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8422T-100	60mL	ST 16	PTFE 2mm	1	\$112.00
LG-8422T-102	125mL	ST 22	PTFE 2mm	1	\$124.76
LG-8422T-104	250mL	ST 22	PTFE 4mm	1	\$147.21
LG-8422T-106	500mL	ST 27	PTFE 4mm	1	\$169.54
LG-8422T-108	1000mL	ST 27	PTFE 4mm	1	\$225.59
LG-8422T-110	2000mL	ST 38	PTFE 6mm	1	\$308.55
LG-8422T-112	4000mL	ST 38	PTFE 8mm	1	\$425.31



Squibb Separatory Funnels with Glass Plug/Stopper and Drip Tip

These Squibb separatory funnels come with glass stopcock plug and glass pennyhead stopper. The bottom of the funnel has a full length Standard Taper inner joint with a drip tip as the delivery system.

CATALOG NO.	CAPACITY	JOINT	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8430-100	60mL	ST 24/40	ST 16	ST 2mm	1	\$169.77
LG-8430-102	125mL	ST 24/40	ST 22	ST 2mm	1	\$169.93
LG-8430-104	250mL	ST 24/40	ST 22	ST 4mm	1	\$190.32
LG-8430-106	500mL	ST 24/40	ST 27	ST 4mm	1	\$225.59
LG-8430-108	1000mL	ST 29/42	ST 27	ST 4mm	1	\$355.54
LG-8430-110	2000mL	ST 29/42	ST 38	ST 6mm	1	\$441.14



Globe Style Separation Funnels with Glass Stopper and Stopcock

This borosilicate glass pear-shaped separatory funnel is supplied with a Standard Taper glass stopper and a Standard Taper glass stopcock.

CATALOG NO.	CAPACITY	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8420-100	30ml	ST 13	ST 2	1	\$84.11
LG-8420-102	60mL	ST 16	ST 2	1	\$85.11
LG-8420-104	125mL	ST 22	ST 2	1	\$99.16
LG-8420-106	250mL	ST 22	ST 4	1	\$112.96
LG-8420-108	500mL	ST 27	ST 4	1	\$138.83
LG-8420-110	1000mL	ST 27	ST 4	1	\$204.35
LG-8420-112	2000mL	ST 38	ST 6	1	\$273.85
LG-8420-114	3000mL	ST 38	ST 8	1	\$350.42
LG-8420-116	4000mL	ST 38	ST 8	1	\$418.68
LG-8420-118	6000ml	ST 38	ST 8	1	\$489.35



Squibb Separatory Funnel with PTFE Plug

This borosilicate glass pear-shaped separatory funnel is supplied as a funnel only or with Standard Taper glass stopper/stopcock and a PTFE plug.

CATALOG NO.	DESCRIPTION	CAPACITY	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8421T-100	Complete	30mL	ST 13	PTFE 2	1	\$96.92
LG-8421T-102	Complete	60mL	ST 16	PTFE 2	1	\$101.80
LG-8421T-104	Complete	125mL	ST 22	PTFE 2	1	\$110.81
LG-8421T-106	Complete	250mL	ST 22	PTFE 4	1	\$132.12
LG-8421T-108	Complete	500mL	ST 27	PTFE 4	1	\$154.61
LG-8421T-110	Complete	1000mL	ST 27	PTFE 4	1	\$206.15
LG-8421T-112	Complete	2000mL	ST 38	PTFE 6	1	\$299.49
LG-8421T-114	Complete	3000mL	ST 38	PTFE 8	1	\$365.73
LG-8421T-116	Complete	4000mL	ST 38	PTFE 8	1	\$432.10
LG-8421T-118	Complete	6000mL	ST 38	PTFE 8	1	\$495.92
LG-8421T-200	Funnel Only	30mL	Fits ST 13	Fits PTFE 2	1	\$83.50
LG-8421T-202	Funnel Only	60mL	Fits ST 16	Fits PTFE 2	1	\$85.19
LG-8421T-204	Funnel Only	125mL	Fits ST 22	Fits PTFE 2	1	\$105.08
LG-8421T-206	Funnel Only	250mL	Fits ST 22	Fits PTFE 4	1	\$116.03
LG-8421T-208	Funnel Only	500mL	Fits ST 27	Fits PTFE 4	1	\$145.49
LG-8421T-210	Funnel Only	1000mL	Fits ST 27	Fits PTFE 4	1	\$191.80
LG-8421T-212	Funnel Only	2000mL	Fits ST 38	Fits PTFE 6	1	\$261.95
LG-8421T-214	Funnel Only	3000mL	Fits ST 38	Fits PTFE 8	1	\$310.60
LG-8421T-216	Funnel Only	4000mL	Fits ST 38	Fits PTFE 8	1	\$366.97
LG-8421T-218	Funnel Only	6000mL	Fits ST 38	Fits PTFE 8	1	\$421.13
ML-1441T-700	Complete	10mL	ST 14/20	PTFE 2	1	\$103.25
ML-1441T-702	Complete	15mL	ST 14/20	PTFE 2	1	\$103.25
ML-1441T-704	Complete	25mL	ST 14/20	PTFE 2	1	\$105.19
ML-1441T-706	Complete	30mL	ST 14/20	PTFE 2	1	\$105.19
ML-1441T-708	Complete	50mL	ST 14/20	PTFE 2	1	\$103.25



Graduated Squibb Cylindrical Separatory Funnel with Glass Stopper/PTFE Plug

The volume of these Squibb separatory funnels is calibrated to the stopcock and offered with increasing/decreasing volume graduations. They feature a Standard Taper inner joint and drip tip at bottom and are supplied with Standard Taper glass stopper and stopcock with PTFE plug.

CATALOG NO.	CAPACITY	SUBDIVISIONS	LOWER JOINT	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8471T-100	60mL	1mL	ST 24/40	ST 16	PTFE 2mm	1	\$234.86
LG-8471T-102	125mL	1mL	ST 24/40	ST 22	PTFE 2mm	1	\$237.01
LG-8471T-104	250mL	5mL	ST 24/40	ST 22	PTFE 2mm	1	\$253.01
LG-8471T-106	250mL	5mL	ST 29/42	ST 22	PTFE 2mm	1	\$293.57
LG-8471T-108	500mL	5mL	ST 24/40	ST 27	PTFE 4mm	1	\$306.68



Globe-Style Separatory Funnel with Glass Plug/Glass Stopper

Globe-shaped separatory funnels have a long stem at bottom and Standard Taper outer neck at top. Supplied with Standard Taper glass stopcock plug and glass pennyhead stopper.


CATALOG NO.	CAPACITY	STEM LENGTH	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8370-100	30mL	175mm	ST 16	ST 2mm	1	\$86.94
LG-8370-102	60mL	175mm	ST 16	ST 2mm	1	\$91.67
LG-8370-104	125mL	175mm	ST 22	ST 4mm	1	\$110.14
LG-8370-108	500mL	175mm	ST 27	ST 4mm	1	\$182.23
LG-8370-110	1000mL	175mm	ST 27	ST 4mm	1	\$247.22
LG-8370-112	2000mL	135mm	ST 38	ST 6mm	1	\$397.46



Globe-Style Separatory Funnel with PTFE Plug/Glass Stopper

Globe-shaped separatory funnels have a long stem at bottom and Standard Taper outer neck at top. Supplied with PTFE stopcock plug and glass pennyhead stopper.


CATALOG NO.	CAPACITY	STEM LENGTH	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8371T-100	30mL	175mm	ST 16	PTFE 2mm	1	\$91.44
LG-8371T-102	60mL	175mm	ST 16	PTFE 2mm	1	\$102.62
LG-8371T-104	125mL	175mm	ST 22	PTFE 2mm	1	\$120.91
LG-8371T-106	250mL	175mm	ST 22	PTFE 2mm	1	\$142.81
LG-8371T-108	500mL	175mm	ST 27	PTFE 4mm	1	\$165.41
LG-8371T-110	1000mL	175mm	ST 27	PTFE 4mm	1	\$237.02
LG-8371T-112	2000mL	130mm	ST 38	PTFE 6mm	1	\$341.68
LG-8371T-116	6000mL	135mm	ST 38	PTFE 8mm	1	\$517.56



Graduated Squibb Cylindrical Separatory Funnel with Glass Plug/Stopper

The volume of these Squibb separatory funnels is calibrated to the stopcock and offered with increasing/decreasing volume graduations. They feature a Standard Taper inner joint and drip tip at bottom and are supplied with Standard Taper glass stopper and glass stopcock.

CATALOG NO.	CAPACITY	SUBDIVISIONS	LOWER JOINT	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8480-100	125mL	1mL	ST 24/40	ST 22	ST 2mm	1	\$232.03
LG-8480-102	250mL	5mL	ST 24/40	ST 22	ST 2mm	1	\$233.06
LG-8480-106	500mL	5mL	ST 24/40	ST 27	ST 4mm	1	\$267.37



Graduated Cylindrical Separatory Funnels with PTFE Plug/Glass Stopper

The volume of these separatory funnels is calibrated to the stopcock and offered with increasing/decreasing volume graduations. Supplied with PTFE plug and Standard Taper glass stopper.

CATALOG NO.	CAPACITY	SUBDIVISIONS	STOPPER SIZE	STOPCOCK SIZE	QUANTITY / PACK	PRICE
LG-8441T-100	125mL	1mL	ST 22	PTFE 2mm	1	\$196.62
LG-8441T-102	250mL	5mL	ST 22	PTFE 2mm	1	\$211.36
LG-8441T-104	500mL	5mL	ST 27	PTFE 4mm	1	\$266.56
LG-8441T-106	1000mL	10mL	ST 27	PTFE 4mm	1	\$267.50
LG-8441T-108	2000mL	20mL	ST 27	PTFE 6mm	1	\$418.91
LG-8441T-200	125mL	1mL	ST 24/40	PTFE 2mm	1	\$235.33
LG-8441T-202	250mL	5mL	ST 24/40	PTFE 2mm	1	\$241.34
LG-8441T-204	500mL	5mL	ST 24/40	PTFE 4mm	1	\$309.26
LG-8441T-206	1000mL	10mL	ST 24/40	PTFE 4mm	1	\$317.49
LG-8441T-208	2000mL	20mL	ST 24/40	PTFE 6mm	1	\$430.63



Standard Taper Pennyhead Stoppers

Product is used extensively with volumetric flasks and separator funnels.

CATALOG NO.	DESCRIPTION	STOPPER SIZE	QUANTITY / PACK	PRICE
LG-10320-104	Glass Pennyhead Stopper	16	1	\$10.47
LG-10320-108	Glass Pennyhead Stopper	22	1	\$13.24
LG-10320-110	Glass Pennyhead Stopper	27	1	\$17.15
LG-10320-114	Glass Pennyhead Stopper	38	1	\$32.09

WEIGHING FUNNELS

Weighing Funnel

Made of borosilicate glass, this funnel provides a means of direct transfer to volumetric flasks, beakers, etc. The base of the funnel is flat to provide stability on the balance pan.



CATALOG NO.	LENGTH	O.D.	QUANTITY / PACK	PRICE
LG-8315-100	40mm	15mm	1	\$20.79
LG-8315-102	65mm	25mm	1	\$24.86