

## VOLUMETRIC CLASS A PIPETTES



### Volumetric Color Coded Class A Pipettes

Calibrated "to deliver" (TD) at 20°C, these borosilicate glass, color coded pipettes are fabricated to meet the requirements for accuracy, appearance, glass quality and calibration line accuracy.

Reference: Volumetric, Class A ASTM E-969

\*NOTE: LG-9350-136 cannot be manufactured to meet Class A specs

CATALOG NO.	CAPACITY	VOLUME TOLERANCE	COLOR	QUANTITY / PACK	PRICE
LG-9350-100	40mL	±0.05mL	White Band	1	\$18.76
LG-9350-101	25mL	±0.03mL	Blue Band	1	\$20.97
LG-9350-102	100mL	±0.08mL	Yellow Band	1	\$20.97
LG-9350-104	9mL	±0.02mL	Black Band	1	\$21.37
LG-9350-106	30mL	±0.03mL	Black Band	1	\$21.37
LG-9350-108	20mL	±0.03mL	Yellow Band	1	\$21.00
LG-9350-110	15mL	±0.03mL	Green Band	1	\$22.40
LG-9350-112	10mL	±0.02mL	Red Band	1	\$21.85
LG-9350-114	80mL	±0.08mL	Blue Band	1	\$22.40
LG-9350-116	6mL	±0.01mL	2 Orange Bands	1	\$22.40
LG-9350-118	8mL	±0.02mL	Blue Band	1	\$24.36
LG-9350-119	7mL	±0.01mL	2 Green Bands	1	\$34.40
LG-9350-120	5mL	±0.01mL	White Band	1	\$26.59
LG-9350-122	1mL	±0.006mL	Blue Band	1	\$26.81
LG-9350-124	2mL	±0.006mL	Orange Band	1	\$28.64
LG-9350-126	4mL	±0.01mL	2 Red Bands	1	\$37.11
LG-9350-128	200mL	Cannot be mfg to meet Class A specs	Blue Band	1	\$37.97
LG-9350-134	0.5mL	±0.006mL	2 Black Bands	1	\$51.50
LG-9350-136	12.5mL	±0.03mL	Orange Band	1	\$93.06

## VOLUMETRIC LAB GRADE PIPETTES



### Pipette-Type Tilting Dispensers

These fixed volume, pipette-type dispensers are designed to deliver measured volumes quickly with reproducibility. To fill, tilt backward so head is vertical. When assembly is returned to an upright position, excess liquid in the head flows into the flask. Tilting forward delivers the measured volume. Complete assembly is supplied with joint clamp (LG-1004) and flask (LG-7750).

CATALOG NO.	DESCRIPTION	CAPACITY	FLASK CAPACITY	JOINT	QUANTITY / PACK	PRICE
LG-7915-108	Tilting Dispenser	5mL	250mL	ST 24/40	1	\$203.67
LG-7915-110	Tilting Dispenser	10mL	500mL	ST 24/40	1	\$205.66
LG-7915-112	Tilting Dispenser	15mL	500mL	ST 24/40	1	\$208.96
LG-7915-118	Tilting Dispenser	25mL	1000mL	ST 24/40	1	\$224.91
LG-7915-132	Tilting Dispenser	50mL	2000mL	ST 24/40	1	\$229.33
LG-7915-138	Tilting Dispenser	100mL	2000mL	ST 34/45	1	\$275.92
LG-7915-208	Tilting Dispenser Head	5mL		ST 24/40	1	\$170.47
LG-7915-210	Tilting Dispenser Head	10mL		ST 24/40	1	\$170.47
LG-7915-212	Tilting Dispenser Head	15mL		ST 24/40	1	\$174.03
LG-7915-218	Tilting Dispenser Head	25mL		ST 24/40	1	\$178.42
LG-7915-232	Tilting Dispenser Head	50mL		ST 24/40	1	\$188.47
LG-7915-238	Tilting Dispenser Head	100mL		ST 34/45	1	\$210.40



### Mohr Class B Measuring Pipettes with White Scales

Calibrated "to deliver" at 20°C, these borosilicate glass pipettes have a permanently colored white scale with the zero line at the top. Capacity of unit is just above the tapered delivery tip. Diameter is reduced at the top for ease in manipulation.

Reference: ASTM E 1293 Style, Type 1, Class B

NOTE: Do not pipette by mouth. Please use a mechanical device.

CATALOG NO.	CAPACITY	SUBDIVISIONS	VOLUME TOLERANCE	COLOR	QUANTITY / PACK	PRICE
LG-9287-110	2mL	1/100mL	±0.02mL	2-White Bands	1	\$10.78
LG-9287-112	5mL	1/10mL	±0.04mL	Blue Band	1	\$10.60



## Serological Color Coded Class B Pipettes

Calibrated "to deliver" (TD) at 20°C, these borosilicate glass pipettes have a permanently colored amber scale with the zero line at the top. The small amount remaining in the tip after completion of free delivery is blown out and added to the main delivery volume (ground band indicator at top). The 5, 10, and 25mL pipettes have openings of reduced diameter.  
Reference: ASTM E 1044, Lab Grade

NOTE: Do not pipette by mouth. A mechanical device is suggested.

CATALOG NO.	CAPACITY	SUBDIVISIONS	VOLUME TOLERANCE	COLOR	QUANTITY / PACK	PRICE
LG-9335-110	1mL	1/100mL	±0.02mL	Yellow Band	1	\$10.81
LG-9335-112	2mL	1/10mL	±0.02mL	Green Band	1	\$11.02
LG-9335-116	5mL	1/10mL	±0.04mL	Blue Band	1	\$11.67
LG-9335-118	25mL	1/10mL	±0.10mL	White Band	1	\$13.29
LG-9335-120	10mL	1/10mL	±0.06mL	Orange Band	1	\$21.71



## Large Orifice Volumetric Pipettes

Calibrated "to deliver" at 20°C, these borosilicate glass, color coded pipettes have a large orifice and are fabricated to meet the requirements for accuracy, appearance, glass quality, and calibration line accuracy.  
Reference: Volumetric, Laboratory Grade, ASTM E-969

CATALOG NO.	CAPACITY	VOLUME TOLERANCE	COLOR	QUANTITY / PACK	PRICE
LG-9355-100	0.5mL	±0.012mL	2-Black Bands	1	\$16.39
LG-9355-102	1mL	±0.012mL	Blue Band	1	\$16.39
LG-9355-104	2mL	±0.012mL	Orange Band	1	\$16.39
LG-9355-106	3mL	±0.02mL	Black Band	1	\$17.51
LG-9355-110	5mL	±0.02mL	White Band	1	\$17.09
LG-9355-120	10mL	±0.04mL	Red Band	1	\$22.37
LG-9355-122	15mL	±0.06mL	Green Band	1	\$24.83
LG-9355-124	20mL	±0.06mL	Yellow Band	1	\$25.10
LG-9355-126	25mL	±0.06mL	Blue Band	1	\$25.71
LG-9355-132	50mL	±0.10mL	Red Band	1	\$33.47
LG-9355-136	100mL	±0.16mL	Yellow Band	1	\$44.87

## DISPOSABLE PIPETTES

### Long-Tip Sample Transfer NMR Pipettes

Wilmad Long-Tip Sample Transfer NMR Pipettes are designed for easy transfer of liquid samples contained in 5mm OD NMR tubes, long neck volumetric flasks or chromatography columns. Manufactured from high quality ASTM Type 1 Class A glass. Resistant to most organic solvents. Transparency provides easy control of sample loading. Easily attach with a latex bulb (804), sold separately. Manufactured in a clean room. Free from organic and inorganic contamination. Special shaping process ensures a smooth surface to minimize sample loss.



CATALOG NO.	DESCRIPTION	LENGTH	QUANTITY / PACK	PRICE
803A	Long tube NMR pipette	13.75in	100	\$85.79
804	Latex bulb		50	\$16.44